




CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

TITANIUM MINI PLATING SYSTEM



SELF FITTING SCREW

MANUFACTURED FROM HIGH QUALITY TITANIUM AND AVAILABLE IN CROSS HEAD SQUARE HEAD OR HEXAGONAL HEAD

SELF FITTING REALIZES ON EASIER AND FASTER

HANDLING AS THE USE OF THE SCREWDRIVER WITH A SCREW HOLDING DEVICE IS NOT NECESSARY
SELF FITTING REQUIRES A SCREW DRIVER BLADE WITH SELF FITTING FEATURES
SELF FITTING CAN AS WELL BE USED WITH COMMON SCREWDRIVER WITH OR WITH OUT SCREW HOLDING DEVICE, BUT THE SELF FITTING CHARACTER OF THE SCREW WILL THEN BE LOST.
ALL SCREWS IN THE SYSTEM ARE SELF TAPPING SCREWS



SELF DRILLING SCREW

MANUFACTURED FROM HIGH QUALITY TITANIUM ALLOY HAVE FOLLOWING CHARACTERISTICS

SELF DRILLING DOES NOT REQUIRE PREDRILLING AND THEREFORE REDUCES OPERATION TIME.

SELF DRILLING WITH **SELF FITTING** SCREW HEADS.

SELF DRILLING WITH A SHARP SCREW TIP FOR INSERTION.

SELF DRILLING IS SILVER COLOR FOR EASY DIFFERENTIATION.

ALL SCREWS IN SYSTEM ARE **SELF TAPPING**

FITTING DRILLS FOR EACH SYSTEM

| DRILL | SCREW |
|-------------------|-----------------------------|
| 1.0MM \emptyset | 1.2MM \emptyset |
| 1.1MM \emptyset | 1.4MM \emptyset EMERGENCY |
| 1.4MM \emptyset | 1.5MM \emptyset |
| | 1.7MM \emptyset |
| 1.5MM \emptyset | 1.9MM EMERGENCY |
| 1.9MM \emptyset | 2.0MM \emptyset |
| | 2.3MM \emptyset |
| 2.0MM \emptyset | 2.5MM \emptyset EMERGENCY |
| | 2.7MM \emptyset |
| | 3.0MM EMERGENCY |



TITANIUM SCREWS
SELF TAPPING



STANDARD



CROSS HEAD
SELF FITING



CROSS HEAD
SELF FITING/DRILLING



HEXAGONAL
SELF FITING



SQUARE HEAD
SELF FITING

| Ø×MM | CAT NO | CAT NO | CAT NO | CAT NO | CAT NO |
|----------|-----------|-----------|-----------|-----------|-----------|
| 2.0×4.0 | CMF-2004A | CMF-2004B | CMF-2004C | CMF-2004D | CMF-2004E |
| 2.0×5.0 | CMF-2005A | CMF-2005B | CMF-2005C | CMF-2005D | CMF-2005E |
| 2.0×6.0 | CMF-2006A | CMF-2006B | CMF-2006C | CMF-2006D | CMF-2006E |
| 2.0×7.0 | CMF-2007A | CMF-2007B | CMF-2007C | CMF-2007D | CMF-2007E |
| 2.0×8.0 | CMF-2008A | CMF-2008B | CMF-2008C | CMF-2008D | CMF-2008E |
| 2.0×9.0 | CMF-2009A | CMF-2009B | CMF-2009C | CMF-2009D | CMF-2009E |
| 2.0×10.0 | CMF-2010A | CMF-2010B | CMF-2010C | CMF-2010D | CMF-2010E |
| 2.0×11.0 | CMF-2011A | CMF-2011B | CMF-2011C | CMF-2011D | CMF-2011E |
| 2.0×12.0 | CMF-2012A | CMF-2012B | CMF-2012C | CMF-2012D | CMF-2012E |
| 2.0×13.0 | CMF-2013A | CMF-2013B | CMF-2013C | CMF-2013D | CMF-2013E |
| 2.0×14.0 | CMF-2014A | CMF-2014B | CMF-2014C | CMF-2014D | CMF-2014E |
| 2.0×15.0 | CMF-2015A | CMF-2015B | CMF-2015C | CMF-2015D | CMF-2015E |
| 2.0×16.0 | CMF-2016A | CMF-2016B | CMF-2016C | CMF-2016D | CMF-2016E |
| 2.0×17.0 | CMF-2017A | CMF-2017B | CMF-2017C | CMF-2017D | CMF-2017E |
| 2.0×18.0 | CMF-2018A | CMF-2018B | CMF-2018C | CMF-2018D | CMF-2018E |
| 2.0×19.0 | CMF-2019A | CMF-2019B | CMF-2019C | CMF-2019D | CMF-2019E |
| 2.0×20.0 | CMF-2020A | CMF-2020B | CMF-2020C | CMF-2020D | CMF-2020E |
| 2.0×21.0 | CMF-2021A | CMF-2021B | CMF-2021C | CMF-2021D | CMF-2021E |
| 2.0×22.0 | CMF-2022A | CMF-2022B | CMF-2022C | CMF-2022D | CMF-2022E |
| 2.0×23.0 | CMF-2023A | CMF-2023B | CMF-2023C | CMF-2023D | CMF-2023E |
| 2.0×24.0 | CMF-2024A | CMF-2024B | CMF-2024C | CMF-2024D | CMF-2024E |



TITANIUM EMERGENCY SCREWS
SELF TAPPING



STANDARD



CROSS HEAD



HEXAGONAL HEAD



HEXAGONAL HEAD

SELF FITTING

SELF FITTING

SELF FITTING

Ø×MM

CAT NO

CAT NO

CAT NO

CAT NO

2.3MM×4.0MM
 2.3MM×5.0MM
 2.3MM×6.0MM
 2.3MM×7.0MM
 2.3MM×8.0MM
 2.3MM×9.0MM
 2.3MM×10.0MM
 2.3MM×11.0MM
 2.3MM×12.0MM
 2.3MM×13.0MM
 2.3MM×14.0MM
 2.3MM×15.0MM
 2.3MM×16.0MM
 2.3MM×17.0MM
 2.3MM×18.0MM
 2.3MM×19.0MM
 2.3MM×20.0MM
 2.3MM×21.0MM
 2.3MM×22.0MM
 2.3MM×23.0MM
 2.3MM×24.0MM

CMF-2304A
 CMF-2305A
 CMF-2306A
 CMF-2307A
 CMF-2308A
 CMF-2309A
 CMF-2310A
 CMF-2311A
 CMF-2312A
 CMF-2313A
 CMF-2314A
 CMF-2315A
 CMF-2316A
 CMF-2317A
 CMF-2318A
 CMF-2319A
 CMF-2320A
 CMF-2321A
 CMF-2322A
 CMF-2323A
 CMF-2324A

CMF-2304B
 CMF-2305B
 CMF-2306B
 CMF-2307B
 CMF-2308B
 CMF-2309B
 CMF-2310B
 CMF-2311B
 CMF-2312B
 CMF-2313B
 CMF-2314B
 CMF-2315B
 CMF-2316B
 CMF-2317B
 CMF-2318B
 CMF-2319B
 CMF-2320B
 CMF-2321B
 CMF-2322B
 CMF-2323B
 CMF-2324B


CMF-2304C
 CMF-2305C
 CMF-2306C
 CMF-2307C
 CMF-2308C
 CMF-2309C
 CMF-2310C
 CMF-2311C
 CMF-2312C
 CMF-2313C
 CMF-2314C
 CMF-2315C
 CMF-2316C
 CMF-2317C
 CMF-2318C
 CMF-2319C
 CMF-2320C
 CMF-2321C
 CMF-2322C
 CMF-2323C
 CMF-2324C

CMF-2304D
 CMF-2305D
 CMF-2306D
 CMF-2307D
 CMF-2308D
 CMF-2309D
 CMF-2310D
 CMF-2311D
 CMF-2312D
 CMF-2313D
 CMF-2314D
 CMF-2315D
 CMF-2316D
 CMF-2317D
 CMF-2318D
 CMF-2319D
 CMF-2320D
 CMF-2321D
 CMF-2322D
 CMF-2323D
 CMF-2324D



PAKISTAN

MINI PLATES



STRAIGHT PLATES

| | | | | | | | |
|-----------|-------------------------|----------------------------|--------------------------|-------------------------|--------------------------|-------------------------|------------------------|
| 2 HOLES | 4 HOLES | 4 HOLES | 4 HOLES | 6 HOLES | 6 HOLES | 8 HOLES | 16 HOLES |
| CMF-10102 | W/O BRIDGE CMF-10104 | MEDIUM BRIDGE CMF-10184 | WITH BRIDGE CMF-10194 | W/O BRIDGE CMF-10106 | WITH BRIDGE CMF-10196 | W/O BRIDGE CMF-10108 | W/O BRIDGE CMF-1016 |

L-PLATES

| | | | | | | | |
|----------------------------------|---------------------------------|------------------------------------|-----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|----------------------------------|
| RIGHT W/O BRIDGE CMF-10104 | LEFT W/O BRIDGE CMF-10224 | RIGHT SHORT BRIDGE CMF-10254 | LEFT SHORT BRIDGE CMF-10264 | RIGHT WITH BRIDGE CMF-10274 | LEFT WITH BRIDGE CMF-10294 | RIGHT LONG BRIDGE CMF-10234 | LEFT LONG BRIDGE CMF-10244 |
|----------------------------------|---------------------------------|------------------------------------|-----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|----------------------------------|



Y-PLATES

DOUBLE Y-PLATES

NOSE PLATES



W/O BRIDGE
CMF-10705



WITH BRIDGE
CMF-10795



W/O BRIDGE
CMF-10606



WITH HOLE BRIDGE
CMF-10607



WITH BRIDGE
CMF-10696



LONG BRIDGE
CMF-10698



Y-PLATE
CMF-10910



20 HOLES
W/O BRIDGE
CMF-10120



40 HOLES
W/O BRIDGE
CMF-10140

ORBITA PLATES



SEMICIRCULAR SHORT
CMF-10806



SEMICIRCULAR LONG
CMF-10808

T-PLATES



90°
W/O BRIDGE
CMF-10306



90°
WITH BRIDGE
CMF-10375



80° RIGHT
WITH BRIDGE
CMF-10355



80° LEFT
WITH BRIDGE
CMF-10365



90°
W/O BRIDGE
CMF-10316



90°
WITH BRIDGE
CMF-10314



L-PLATES



5 HOLES
RIGHT LONG BRIDGE
CMF-10275



5 HOLES
LEFT LONG BRIDGE
CMF-10295



LEFT
112°
CMF-10210



RIGHT
112°
CMF-10211



RIGHT WITH
BRIDGE 112°
CMF-10212



LEFT WITH
BRIDGE 112°
CMF-10213


QUADRATIC/RECTANGULAR PLATES



H-PLATE
CMF-10912



CMF-10504



CMF-10506



CMF-10404



CMF-10406



Ortho Innovations.

PAKISTAN

INSTRUMENTATION

DRILL GUIDE WITH HANDLE
2.0MM SCREWS



WITH INSERT SINGLE HANDLE
CMF-20318 CMF-20218

DEPTH GAUGE
19CM, 7"
CMF-20018




SPIRAL DRILL BIT
1.5MMØ
7.0MM STOP×50MM CMF-15050
7.0MM STOP×70MM CMF-15070
22MM STOP ×50MM CMF-15950
22MM STOP ×70MM CMF-15970
22MM STOP × 105MM CMF-15105



TRANSBUCCALE DRILL GUIDE

COMPLETE FOR SCREW DRIVER W/O SCREW HOLDING DEVICE
3.0MMØ CMF-30000
HANDLE WITH SLEVE CMF-30010
CHEK-HOLDING CMF-30020
TROCER CMF-30030
DRILL GUIDE CMF-30040



COMPLETE FOR SCREW DRIVER WITH SCREW HOLDING DEVICE
6.0MMØ CMF-6000
HANDLE WITH SLEVE CMF-60010
CHEEK – HOLDING CMF-60020
TROCER CMF-60030
DRILLGUIDE CMF-60040



Ortho Innovations.

PAKISTAN





PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7" CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313



Ortho Innovations.

PAKISTAN

PLATE HOLDING CLAMP
14CM, 5-1/2" CMF-50013



PLATE HOLDING FORK
18CM, 7" CMF-2218



PLATE HOLDING FORCEPS STRAIGHT
18CM, 7" CMF-40018



PLATE HOLDING FORCEPS CURVED
18CM, 7" CMF-45018





Orthon Innovations.

PAKISTAN



| | | | | |
|---|--|--|--|---|
| SCREW DRIVER DESMANTABLE COMPLETE HANDLE SCREW HOLDING DEVICE BLADE | STANDARD CMF-20118 CMF-00000 CMF-20000 CMF-20100 | CROSS HEAD CMF-20218 CMF-00000 CMF-20000 CMF-20200 | HEXAGONAL HEAD CMF-20318 CMF-00000 CMF-20000 CMF-20300 | SQUARE HEAD CMF-20418 CMF-00000 CMF-20000 CMF-20400 |
|---|--|--|--|---|



| | | | | |
|--|----------|--------------------------|------------------------------|---------------------------|
| SCREWDRIVER DESMENTABLE COMPLETE SELF FITTING | STANDARD | CROSS HEAD CMF-202185 | HEXAGONAL HEAD CMF-203185 | SQUARE HEAD CMF-204185 |
| HANDLE BLADE | | CMF-00000 CMF-202005 | CMF-00000 CMF-203005 | CMF-00000 CMF-204005 |




| | | | | |
|------------------------------|-----------------------|-------------------------|-----------------------------|--------------------------|
| SCREW DRIVER COMPLETE | STANDARD CMF-20118 | CROSS HEAD CMF-20218 | HEXAGONAL HEAD CMF-20318 | SQUARE HEAD CMF-20418 |
| SCREW DRIVER SELF FITTING | | CMF-202185 | CMF-203185 | CMF-205185 |



CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

THE LOW PROFILE SYSTEM

TITANIUM SCREWS SELF-TAPPING



Ø×MM

STANDARD


CROSS HEAD (SEFL FITTING)

- 1.5×4.0
- 1.5×5.0
- 1.5×6.0
- 1.5×7.0
- 1.5×8.0
- 1.5×9.0
- 1.5×10.0
- 1.5×11.0
- 1.5×12.0
- 1.5×13.0
- 1.5×14.0
- 1.5×15.0
- 2.0×4.0
- 2.0×5.0
- 2.0×6.0
- 2.0×7.0
- 2.0×8.0
- 2.0×9.0
- 2.0×10.0
- 2.0×11.0
- 2.0×12.0
- 2.0×13.0
- 2.0×14.0
- 2.0×15.0
- 2.0×16.0
- 2.0×17.0
- 2.0×18.0
- 2.0×19.0
- 2.0×20.0
- 2.0×21.0
- 2.0×22.0
- 2.0×23.0
- 2.0×24.0

- CMF-15104
- CMF-15105
- CMF-15106
- CMF-15107
- CMF-15108
- CMF-15109
- CMF-15110
- CMF-15111
- CMF-15112
- CMF-15113
- CMF-15114
- CMF-15115
- CMF-20704
- CMF-20705
- CMF-20706
- CMF-20707
- CMF-20708
- CMF-20709
- CMF-20710
- CMF-20711
- CMF-20712
- CMF-20713
- CMF-20714
- CMF-20715
- CMF-20716
- CMF-20717
- CMF-20718
- CMF-20719

- CMF-15204
- CMF-15205
- CMF-15206
- CMF-15207
- CMF-15208
- CMF-15209
- CMF-15210
- CMF-15211
- CMF-15212
- CMF-15213
- CMF-15214
- CMF-15215
- CMF-20204
- CMF-20205
- CMF-20206
- CMF-20207
- CMF-20208
- CMF-20209
- CMF-20210
- CMF-20211
- CMF-20212
- CMF-20213
- CMF-20214
- CMF-20215
- CMF-20216
- CMF-20218
- CMF-20218
- CMF-20219
- CMF20220
- CMF-20221
- CMF-20222
- CMF-20223
- CMF-20224


LOW PROFILE SYSTEM



HEIGHT 1.1MM

THICKNESS 0.6MM

MINI PLATING SYSTEM




HEIGHT 1.3MM

THICKNESS 1.0MM



LOW PROFILE PLATES

STRAIGHT PLATES



| | | | | | | |
|----------------------|------------------------------------|--|-------------------------------------|------------------------------------|-------------------------------------|------------------------------------|
| 2 HOLES CMF-06102 | 4 HOLES W/O BRIDGE CMF-06104 | 4 HOLES MEDIUM BRIDGE CM F-06184 | 4 HOLES WITH BRIDGE CMF-06194 | 6 HOLES W/O BRIDGE CMF-06106 | 6 HOLES WITH BRIDGE CMF-06196 | 8 HOLES W/O BRIDGE CMF-06108 |
|----------------------|------------------------------------|--|-------------------------------------|------------------------------------|-------------------------------------|------------------------------------|

L-PLATES



| | | | | | | | |
|----------------------------------|---------------------------------|------------------------------------|-----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|----------------------------------|
| RIGHT W/O BRIDGE CMF-06204 | LEFT W/O BRIDGE CMF-06224 | RIGHT SHORT BRIDGE CMF-06254 | LEFT SHORT BRIDGE CMF-06264 | RIGHT WITH BRIDGE CMF-06274 | LEFT WITH BRIDGE CMF-06294 | RIGHT LONG BRIDGE CMF-06275 | LEFT LONG BRIDGE CMF-06295 |
|----------------------------------|---------------------------------|------------------------------------|-----------------------------------|-----------------------------------|----------------------------------|-----------------------------------|----------------------------------|

T-PLATES

DOUBLE T-PLATES



| | | | | | |
|-------------------------------|---------------------------------|---------------------------------------|--------------------------------------|--------------------------------|---------------------------------|
| 90° W/O BRIDGE CMF-6306 | 90° WITH BRIDGE CMF-06375 | 80° RIGHT WITH BRIDGE CMF-06355 | 80° LEFT WITH BRIDGE CMF-06365 | 90° W/O BRIDGE CMF-06316 | 90° WITH BRIDGE CMF-06314 |
|-------------------------------|---------------------------------|---------------------------------------|--------------------------------------|--------------------------------|---------------------------------|

ORBITA PLATES



| | |
|---------------------------------|--------------------------------|
| SEMICIRCULAR SHORT CMF-06806 | SEMICIRCULAR LONG CMF-06808 |
|---------------------------------|--------------------------------|



PAKISTAN

LOW PROFILE PLATES



Y-PLATES

DOUBLE Y-PLATES



W/O BRIDGE
CMF-06705



WITH BRIDGE
CMF-06795



W/O BRIDGE
CMF-06606



WITH BRIDGE
CMF-06696



WITH LONG BRIDGE
CMF-06698

NOSE PLATES



H-PLATE
CMF-06912



Y-PLATE
CMF-06910



BURR HOLE COVER
Ø 7MM CMF-06990
Ø 10MM CMF-06995
Ø 14MM CMF-06997

QUADRATIC AND RECTANGULAR PLATES



RECTANGULAR
QUADRATIC
CMF-06404
CMF-06504



CMF-06406
CMF-06506



CMF-06408
CMF-06508



CMF-06416
CMF-06416



CMF-06424
CMF-06524




Ortho Innovations.


PAKISTAN



14x2HOLES
CMF-06528



9x14 HOLES
CMF-06599



35x12 HOLES
CMF-06550



Ortho Innovations.

PAKISTAN





PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7" CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313



Ortho Innovations.

PAKISTAN

PLATE HOLDING CLAMP
14CM, 5-1/2" CMF-50013



PLATE HOLDING FORK
18CM, 7" CMF-2218



PLATE HOLDING FORCEPS STRAIGHT
18CM, 7" CMF-40018



PLATE HOLDING FORCEPS CURVED
18CM, 7" CMF-45018





| | | | | |
|--|--|--|--|---|
| SCREW DRIVER DESMANTABE COMPLETE HANDLE SCREW HOLDING DEVICE BLADE | STANDARD CMF-20118 CMF-00000 CMF-20000 CMF-20100 | CROSS HEAD CMF-20218 CMF-00000 CMF-20000 CMF-20200 | HEXAGONAL HEAD CMF-20318 CMF-00000 CMF-20000 CMF-20300 | SQUARE HEAD CMF-20418 CMF-00000 CMF-20000 CMF-20400 |
|--|--|--|--|---|



| | | | | |
|--|----------|--------------------------|------------------------------|---------------------------|
| SCREWDRIVER DESMENTABLE COMPLETE SELF FITTING | STANDARD | CROSS HEAD CMF-202185 | HEXAGONAL HEAD CMF-203185 | SQUARE HEAD CMF-204185 |
| HANDLE BLADE | | CMF-00000 CMF-202005 | CMF-00000 CMF-203005 | CMF-00000 CMF-204005 |



| | | | | |
|------------------------------|-----------------------|-------------------------|-----------------------------|--------------------------|
| SCREW DRIVER COMPLETE | STANDARD CMF-20118 | CROSS HEAD CMF-20218 | HEXAGONAL HEAD CMF-20318 | SQUARE HEAD CMF-20418 |
| SCREW DRIVER SELF FITTING | | CMF-202185 | CMF-203185 | CMF-205185 |



Ortho Innovations

PAKISTAN

CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

MEDIUM PLATING SYSTEM

TITANIUM SCREWS SELF TAPPING



STANDARD



CROSS HEAD
SELF FITING



CROSS HEAD
SELF FITING/DRILLING



SQUARE HEAD
SELF FITING

Ø×MM
1.5×3.0
1.5×4.0
1.5×5.0
1.5×6.0
1.5×7.0
1.5×8.0
1.5×9.0
1.5×10.0
1.5×11.0
1.5×11.0
1.5×12.0
1.5×13.0
1.5×14.0
1.5×15.0
1.5×16.0
1.5×17.0
1.5×18.0
1.5×19.0
1.5×20.0
1.5×21.0
1.5×22.0
1.5×23.0
1.5×24.0

CAT NO
CMF-15103
CMF-15104
CMF-15105
CMF-15106
CMF-15107
CMF-15108
CMF-15109
CMF-15110
CMF-15111
CMF-15112
CMF-15113
CMF-15114
CMF-15115
CMF-15116
CMF-15117
CMF-15118
CMF-15119
CMF-15120

CAT NO
CMF-15203
CMF-15204
CMF-15205
CMF-15206
CMF-15207
CMF-15208
CMF-15209
CMF-15210
CMF-15211
CMF-15212
CMF-15213
CMF-15214
CMF-15215
CMF-15216
CMF-15217
CMF-15218
CMF-15219
CMF-15220
CMF-15221
CMF-15222
CMF-15223
CMF-15224

CAT NO
CMF-15303
CMF-15304
CMF-15305
CMF-15306
CMF-15307
CMF-15308
CMF-15309
CMF-15310
CMF-15311

CAT NO
CMF-15403
CMF-15404
CMF-15405
CMF-15406
CMF-15407
CMF-15408
CMF-15409
CMF-15410
CMF-15411
CMF-15412
CMF-15413
CMF-15414
CMF-15415
CMA-15416
CMF-15417
CMF-15417
CMF-15419
CMF-15420



TITANIUM EMERGENCY SCREWS
SELF TAPPING



STANDARD



CROSS HEAD
SELF FITING



CROSS HEAD
SELF FITING/DRILLING



SQUARE HEAD
SELF FITING

Ø×MM
1.7×3.0
1.7×4.0
1.7×5.0
1.7×6.0
1.7×7.0
1.7×8.0
1.7×9.0
1.7×10.0
1.7×11.0
1.7×12.0
1.7×13.0
1.7×14.0
1.7×15.0
1.7×16.0
1.7×17.0
1.7×18.0
1.7×19.0
1.7×20.0
1.7×21.0
1.7×22.0
1.7×23.0
1.7×24.0
1.7×25.0

CAT NO
CMF-17103
CMF-17104
CMF-17105
CMF-17106
CMF-17107
CMF-17108
CMF-17109
CMF-17110
CMF-17111
CMF-17112
CMF-17113
CMF-17114
CMF-17115
CMF-17116
CMF-17117
CMF-17118
CMF-17119
CMF-17120
CMF-17121
CMF-17122
CMF-17123
CMF-17124
CMF-17125

CAT NO
CMF-17203
CMF-17204
CMF-17205
CMF-17206
CMF-17207
CMF-17208
CMF-17209
CMF-17210
CMF-17211
CMF-17212
CMF-17213
CMF-17214
CMF-17215
CMF-17216
CMF-17217
CMF-17218
CMF-17219
CMF-17220
CMF-17221
CMF-17222
CMF-17223
CMF-17224
CMF-17225

CAT NO
CMF-17303
CMF-17304
CMF-17305
CMF-17306
CMF-17307
CMF-17308
CMF-17309
CMF-17310
CMF-17311

CAT NO
CMF-17403
CMF-17404
CMF-17405
CMF-17406
CMF-17407
CMF-17408
CMF-17409
CMF-17410
CMF-17411
CMF-17412
CMF-17413
CMF-17414
CMF-17415
CMA-17416
CMF-17417
CMF-17419
CMF-17420



TITANIUM EMERGENCY SCREWS
SELF TAPPING



STANDARD



CROSS HEAD
SELF FITING




SQUARE HEAD
SELF FITING

| Ø×MM | CAT NO | CAT NO | CAT NO |
|----------|-----------|-----------|-----------|
| 1.9×4.0 | CMF-19104 | CMF-19204 | CMF-19304 |
| 1.9×5.0 | CMF-19105 | CMF-19205 | CMF-19305 |
| 1.9×6.0 | CMF-19106 | CMF-19206 | CMF-19306 |
| 1.9×7.0 | CMF-19107 | CMF-19207 | CMF-19307 |
| 1.9×8.0 | CMF-19108 | CMF-19208 | CMF-19308 |
| 1.9×9.0 | CMF-19109 | CMF-19209 | CMF-19309 |
| 1.9×10.0 | CMF-19110 | CMF-19210 | CMF-19310 |
| 1.9×11.0 | CMF-19111 | CMF-19211 | CMF-19311 |
| 1.9×12.0 | CMF-19112 | CMF-19212 | CMF-19312 |
| 1.9×13.0 | CMF-19113 | CMF-19213 | CMF-19313 |
| 1.9×14.0 | CMF-19114 | CMF-19214 | CMF-19314 |
| 1.9×15.0 | CMF-19115 | CMF-19215 | CMF-19315 |
| 1.9×16.0 | CMF-19116 | CMF-19216 | CMF-19316 |
| 1.9×17.0 | CMF-19117 | CMF-19217 | CMF-19317 |
| 1.9×18.0 | CMF-19118 | CMF-19218 | CMF-19318 |
| 1.9×19.0 | CMF-19119 | CMF-19219 | CMF-19319 |
| 1.9×20.0 | CMF-19120 | CMF-19220 | CMF-19320 |



MEDIUM PLATES



STRAIGHT PLATES





L-PLATES



RIGHT 90°
W/O BRIDGE
CMF-06204



LEFT 90°
W/O BRIDGE
CMF-06224



RIGHT 90°
SHORT BRIDGE
CMF-06254



LEFT 90°
SHORT BRIDGE
CMF-06264



RIGHT 90°
MED BRIDGE
CMF-06274



LEFT 90°
MEDI BRIDGE
CMF-06294



RIGHT 90°
LARGE BRIDGE
CMF-06234



LEFT 90°
LARGE BRIDGE
CMF-06244



RIGHT
EX LONG BRIDGE
CMF-06275



LEFT
EX LONG BIRDGE
CMF-06295



RIGHT 112°
W/O BRIDGE
CMF-06210



LEFT 112°
W/O BRIDGE
CMF-06211



RIGHT 112°
MEDIUM BRIDGE
CMF-06214



LEFT 112°
MEDIUM BRIDGE
CMF-06215



RIGHT 112°
LONG BRIDGE
CMF-06216



LEFT 112°
LONG BRIDGE
CMF-06217



RIGHT 112°
EXT LONG BRIDGE
CMF-06218



LEFT 112°
EXT LONG BRIDGE
CMF-06219

T-PLATES

DOUBLE- T- PLATE

Y-PLATES



4 HOLES
SHORT BRIDGE
CMF-06304



4 HOLES
MEDIUM BRIDGE
CMF-06305



6 HOLES
W/O BRIDGE
CMF-06306



5 HOLES
WITH BRIDGE
CMF-0606307



6 HOLES
W/O BRIDGE
CMF-06316



5 HOLES
W/O BRIDGE
CMF-06705



5 HOLES
MEDIUM BRIDGE
CMF-06785



5 HOLES
LONG BRIDGE
CMF-06795



DOUBLE- Y- PLATES

ORBITA PLATES



BURR HOLE COVER

ZYGOMA PLATE



QUADRATIC AND RECTANGULAR GRILLE-PLATES





INSTRUMENTATION

SPIRAL DRILL BIT

1.4MMØ

5.0MM STOP×45MM CMF-14005

7.0MM STOP×45MM CMF-14007

9.0MM STOP 45M CMF-14009

15MM STOP ×45MM CMF-14015

SPIRAL DRILL BIT

1.1MM Ø

5MM STOP×45MM CMF-11005

7MM STOP×45MM CMF-11007

9MM STOP×45MM CMF-11009

15MM STOP×45MM CMF-11015



PLATE HOLDING FORCEPS STRAIGHT

18CM, 7" CMF-40018



PLATE HOLDING FORCEPS CURVED

18CM, 7" CMF-45018





Ortho Innovations.

PAKISTAN




PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7" CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313



PLATE HOLDING CLAMP
14CM, 5-1/2" CMF-50013



| | | | |
|----------------------|-----------|------------|-------------|
| SCREW DRIVER | STANDARD | CROSS HEAD | SQUARE HEAD |
| DESMANTABLE COMPLETE | CMF-20117 | CMF-20217 | CMF-20417 |
| HANDLE | CMF-00000 | CMF-00000 | CMF-00000 |
| SCREW HOLDING DEVICE | CMF-20000 | CMF-20000 | CMF-20000 |
| BLADE | CMF-20100 | CMF-20200 | CMF-20400 |



| | | | |
|-------------------------|-----------|------------|-------------|
| SCREWDRIVER DESMENTABLE | STANDARD | CROSS HEAD | SQUARE HEAD |
| COMPLETE SELF FITTING | CMF-20117 | CMF-20217 | CMF-20417 |
| HANDLE | | CMF-00000 | CMF-00000 |
| BLADE | | CMF-20200 | CMF-20400 |



| | | | |
|--------------|-----------|------------|-------------|
| SCREW DRIVER | STANDARD | CROSS HEAD | SQUARE HEAD |
| COMPLETE | CMF-20117 | CMF-20217 | CMF-20417 |
| SCREW DRIVER | | CMF-20217 | CMF-20517 |
| SELF FITTING | | | |



CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

MICRO PLATTING SYSTEM TITANIUM SCREWS SELF-TAPPING



STANDARD



CROSS HEAD
SELF FITTING

TITANIUM EMERGENCY SCREWS SELF-TAPPING



STANDARD



CROSS HEAD
SELF FITTING

Ø×MM

CAT NO

CAT NO

Ø×MM

CAT NO

CAT NO

1.2×2.0
1.2×3.0
1.2×4.0
1.2×5.0
1.2×6.0
1.2×7.0
1.2×8.0
1.2×9.0
1.2×10.0
1.2×11.0
1.2×12.0
1.2×13.0
1.2×14.0
1.2×15.0

CMF-12102
CMF-12103
CMF-12104
CMF-12105
CMF-12106
CMF-12107
CMF-12108
CMF-12109
CMF-12110
CMF-12111
CMF-12112
CMF-12113
CMF-12114
CMF-12115

CMF-12202
CMF-12203
CMF-12204
CMF-12205
CMF-12206
CMF-12207
CMF-12208
CMF-12209
CMF-12210
CMF-12211
CMF-12212
CMF-12213
CMF-12214
CMF-12215


1.4×3.0
1.4×4.0
1.4×5.0
1.4×6.0
1.4×7.0
1.4×8.0
1.4×9.0
1.4×10.0
1.4×11.0
1.4×12.0
1.4×13.0
1.4×14.0
1.4×15.0

CMF-O-14103
CMF-O-14104
CMF-O-14105
CMF-O-14106
CMF-O-14107
CMF-O-14108
CMF-O-14109
CMF-O-14110
CMF-O-14111
CMF-O-14112
CMF-O-14113
CMF-O-14114
CMF-O-14115

CMF-O-14203
CMF-O-14204
CMF-O-14205
CMF-O-14206
CMF-O-14207
CMF-O-14208
CMF-O-14209
CMF-O-14210
CMF-O-14211
CMF-O-14212
CMF-O-14213
CMF-O-14214
CMF-O-14215



MICRO PLATES



PAKISTAN



STRAIGHT PLATES



2 HOLES
CMF-06102



4 HOLES
CMF-06104



6 HOLES
CMF-06106



8 HOLES
CMF-06108



12 HOLES
CMF-06112



16 HOLES
CMF-06116



24 HOLES
CMF-06124

L-PLATES



5 HOLES
RIGHT SHORT
CMF-06205



5 HOLES
LEFT SHORT
CMF-06245



8 HOLES
RIGHT LONG
CMF-0606208



8 HOLES
LEFT LONG
CMF-06248

T-PLATES



90° 5 HOLES
SHORT
CMF-06305



90° 5 HOLES
LONG
CMF-06307

Y-PLATES



6 HOLES
CMF-06706

DOUBLE- Y-PLATES



SHORT 6 HOLES
CMF-06606



LONG 7 HOLES
CMF-06607

H-PLATES



7 HOLES
CMF-06807



9 HOLES
CMF-06809



14 HOLES
CMF-06814



PAKISTAN

QUADRATIC AND RECTANGULAR
GRILLE-PLATES




4 HOLES
CMF-06504



6 HOLES
CMF-06506



8 HOLES
CMF-06508



16 HOLES
CMF-06516



24 HOLES
CMF-06524




4 HOLES
CMF-06404



6 HOLES
CMF-06406



8 HOLES
CMF-06408



16 HOLES
CMF-06416



24 HOLES
CMF-06424

BURR HOLE COVER



7MM Ø
CMF-06997



10MM Ø
CMF-06995



14MM
CMF-06990


ORBITA PLATE



SEMICIRCULAR, 8 HOLES
CMF-06901



TITANIUM MICRO MESH



0.1MM

| | |
|----------|----------|
| CMF-1990 | 90×20MM |
| CMF-1991 | 120×60MM |
| CMF-1992 | 60×60MM |

0.2MM

| | |
|-----------|----------|
| CMF-02990 | 90×120MM |
| CMF-02991 | 120×60MM |
| CMF-02992 | 60×60MM |

SMITH WIRE CUTTING SCISSORS

TC. 16CM, 6"

| | |
|-----------------------|---------------------|
| STRAIGHT CMF-84416 | CURVED CMF-84516 |
|-----------------------|---------------------|

WIRE CUTTING SCISSOR TC

12CM
CMF-83512





INSTRUMENTATION

SPIRAL DRILL BIT

1.0MMØ
4.0MM STOP×46MM CMF-15906
8.0MM STOP×50MM CMF-15907
12.0MM STOP 54M CMF-15908



DEPTH GAUGE

19CM, 7"
CMF-20018



PLATE HOLDING FORCEPS STRAIGHT

18CM, 7" CMF-40018




PLATE HOLDING FORCEPS CURVED

18CM, 7" CMF-45018



PLATE AND WIRE CUTTING PLIER W

12.5CM, 4-3/4" CMF-70412





Ortho Innovations.


PAKISTAN



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313



PLATE HOLDING CLAMP
14CM, 5-1/2" CMF-50013



DESMANTABLEABLE SCREW DRIVER
COMPLETE

| | |
|----------------------------|--------------------------|
| STANDARD CMF-20117 | QCROSS HEAD CMF-20217 |
| HANDLE CMF-00000 | CMF-00000 |
| TENTION PLIER CMF-20000 | CMF-20000 |
| BLADE CMF-20100 | CMF-20200 |



SCREW DRIVER DESMENTABLE
COMPLETE SELF FITTING

| | |
|----------|-------------------------|
| STANDARD | CROSS HEAD CMF-20217 |
| HANDLE | CRF-00000 |
| BLADE | CMF-20200 |



SCREW DRIVER FIXED COMPLETE

| | |
|---------------------------------------|-------------------------|
| STANDARD CMF-20117 | CROSS HEAD CMF-20217 |
| SELF FITTING COMPLETE SCREW DRIVER | CMF-20217 |





Ortho Innovations

PAKISTAN

CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

THE MANDIBULAR RECONSTRUCTION SYSTEM

TITANIUM SCREWS SELF TAPPING



CROSS HEAD

| Ø×MM | CAT NO |
|--------|-----------|
| 2.7×8 | CMF-27298 |
| 2.7×10 | CMF-27210 |
| 2.7×12 | CMF-27212 |
| 2.7×14 | CMF-27214 |
| 2.7×16 | CMF-27216 |
| 2.7×18 | CMF-27218 |
| 2.7×20 | CMF-27220 |
| 2.7×22 | CMF-27222 |

TITANIUM EMERGENCY SCREW SELF TAPPING



CROSS HEAD

| Ø×MM | CAT NO |
|--------|-----------|
| 3.0×10 | CMF-30210 |
| 3.0×12 | CMF-30212 |
| 3.0×14 | CMF-30214 |

TITANIUM FASTENING SCREWS




CROSS HEAD

| CAT NO |
|-----------|
| CMF-28307 |



RECONSTRUCTION PLATES



8 HOLES

RECONSTRUCTION PLATE CMF-22008
TIN BENDING TEMPLATE CMF-22006



11 HOLES

CMF-22011
CMF-22011A



14 HOLES

CMF-22014
CMF-22014A



17 HOLES

CMF-22017
CMF-22017A



4+8 HOLES LEFT
CMF-22312
CMF-22312

RECONSTRUCTION PLATE
TIN BENDING TEMPLATE



4+8 HOLES RIGHT
CMF-22212
CMF-22212


RECONSTRUCTION PLATE
TIN BENDING TEMPLATE



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+11 HOLES LEFT
CMF-22315
CMF-22215

4+11 HOLES RIGHT
CMF-22215
CMF-22215



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+14 HOLES LEFT
CMF-22318
CMF-22218

4+14 HOLES RIGHT
CMF-22218
CMF-22218



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+17+4 HOLES RIGHT

CMF-22425
CMF-22425



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+17+4 HOLES LEFT
CMF-22427
CMF-22427



PROSTHESES



TEMPORARY CONDYLAR PROSTHESES
RIGHT SHORT. CMF-25153




TEMPORARY CONDYLAR PROSTHESES
LEFT SHORT CMF-22150




TEMPORARY CONDYLAR
PROSTHESES LEFT 5 HOLES
CMF-25152



TEMPORARY CONDYLAR
PROSTHESES RIGHT 5 HOLES
CMF-25155



BENDING TEMPLATES
CMF-22000





Ortho Innovations.

PAKISTAN

INSTRUMENTATION

PLATE BENDER 15 CM LENGTHS

CMF-22150



PLATE HOLDING FORCEPS LEFT
15.5CM CMF-25154



PLATE BENDING FORCEPS RIGHT
15.5 CM CMF-25155



PLATE BENDING PLIER
LENGTH 16.5MM FOR
2.2&2.5MM PLATES
CMF-25160





Ortho Innovations.

PAKISTAN



SCREWDRIVER DESMENTABLE
COMPLETE SELF FITTING
CROSS HEAD CMF-27218
HANDLE CMF-00000
BLADE CMF-27200



SCREWDRIVER DESMENTABLE
COMPLETE CMF-27218
HANDLE CMF-00000
BLADE CMF-27200
SCREW HOLDER CMF-27000A



TAP. 2.7MM SCREWS
CMF-27000



DRILL GUIDE (MANDIBULAR SYSTEM)
CMF-28100



SPIRAL DRILL Ø 2.0MM, 35 MM TOP
103MM CMF-20103



Ortho Innovations.

PAKISTAN



DEPTH GAUGE (UP TO 45 MM SCREWS)
CMF-27018



PLATE CUTTING PLIERS 57.5CM LENGTH
CMF-47000

TRANSBUCCALE DRILL GUIDE

| | |
|---|-----------|
| COMPLETE FOR SCREW DRIVER W/O SCREW HOLDING DEVICE | |
| 3.5MMØ | CMF-35000 |
| HANDLE WITH SLEVE | CMF-35010 |
| CHEK-HOLDING | CMF-35020 |
| TROCAR | CMF-35030 |
| DRILL GUIDE | CMF-35040 |
| COMPLETE FOR SCREW DRIVER WITH SCREW HOLDING DEVICE | |
| 6.0MMØ | CMF-60000 |
| HANDLE WITH SLEVE | CMF-60010 |
| CHEEK – HOLDING | CMF-60020 |
| TROCAR | CMF-60030 |
| DRILLGUIDE | CMF-60040 |

