




Ortho Innovations.

PAKISTAN

ORTHON ANGLED BLADE PLATE INSTRUMENTS SET



ABP-001

TRIPPLE DRILL GUIDE 130°



ABP-002

ROUTER



ABP-003

CHIESEL GUIDE





Ortho Innovations

PAKISTAN

ABP-004 CEILING CHISEL



ABP-005 INSERTER/EXTRACTOR



ABP-006 PLATE IMPACTER



ABP-007 SLOTTED HAMMER



ABP-008 POSITIONING PLATE 90/50/40°

ABP-009 POSITIONING PLATE 80/70/30°

ABP-010 POSITIONING PLATE 100/60/20°



ABP-011 ANGLED PLATE GUID





Ortho Innovations.


PAKISTAN

AB-012 NEUTRAL & LOAD DRILL GUIDE



CONDYLAR PLATES 95°

CAT NO	HOLES	BALDE LENGTH	SHAFT LENGTH
ABP-101	5	40MM	92MM
ABP-102	5	50MM	92MM
ABP-103	5	60MM	92MM
ABP-104	5	70MM	92MM
ABP-105	7	40MM	124MM
ABP-106	7	50MM	124MM
ABP-107	7	60MM	124MM
ABP-108	7	70MM	124MM
ABP-109	9	40MM	156MM
ABP-110	9	50MM	156MM
ABP-111	9	60MM	156MM
ABP-112	9	70MM	156MM
ABP-113	12	50MM	204MM
ABP-114	12	60MM	204MM
ABP-115	12	70MM	204MM
ABP-116	12	80MM	204MM
ABP-117	14	50MM	236MM





ABP-118	14	60MM	236MM
ABP-119	14	70MM	236MM
ABP-120	14	80MM	236MM
ABP-121	16	50MM	268MM
ABP-122	16	60MM	268MM
ABP-123	16	70MM	268MM
ABP-124	16	80MM	268MM
ABP-125	18	50MM	300MM
ABP-126	18	60MM	300MM
ABP-127	18	70MM	300MM
ABP-128	18	80MM	300MM

IN TITANIUM THE SAME SPECIFICATION IS FORM
ABP-201 TO ABP-228

PROFILE


THICKNESS	5.7MM
WIDTH	16.2MM
HOLE SPACING	16.0MM
U PROFILE BLADE	6.6/16.2MM

ANGLED BLADE PLATES 130°

ABP-301	4	50MM	60MM
ABP-302	4	60MM	60MM
ABP-303	4	70MM	60MM
ABP-304	4	80MM	60MM
ABP-305	4	85MM	60MM
ABP-306	4	90MM	60MM
ABP-307	4	95MM	60MM
ABP-308	4	100MM	60MM
ABP-309	4	105MM	60MM
ABP-310	4	110MM	60MM
ABP-311	6	50MM	104MM
ABP-312	6	60MM	104MM
ABP-313	6	70MM	104MM
ABP-314	6	80MM	104MM
ABP-315	6	90MM	104MM
ABP-316	9	50MM	152MM
ABP-317	9	60MM	152MM
ABP-318	9	70MM	152MM

PROFILE

THICKNESS	5.4MM
WIDTH	14.0MM
HOLE SPACING	16.0MM
T-PROFILE BLADE	5.0/11.7MM





PAKISTAN

ABP-319	9	80MM	152MM
ABP-320	9	90MM	152MM
ABP-321	12	70MM	200MM
ABP-322	12	80MM	200MM
ABP-323	12	90MM	200MM
ABP-324	14	70MM	232MM
ABP-325	14	80MM	232MM
ABP-326	14	90MM	232MM

IN TITANIUM SAME SPECIFICATION IS FROM ABP-301 TO ABP-326

PROFILE

<i>THICKNESS</i>	<i>5.7MM</i>
<i>WIDTH</i>	<i>16.2MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>
<i>U-PROFILE BLADE</i>	<i>6.6/16.2MM</i>

HIP PLATE 110°

ABP-401	4	65MM
ABP-402	4	75MM
ABP-403	4	85MM

IN TITANIUM SAME SOZES ARE FROM ABP-501 TO ABP-503



PROFILE

<i>THICKNESS</i>	<i>5.7MM</i>
<i>WIDTH</i>	<i>16.2MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>
<i>U-PROFILE BLADE</i>	<i>6.6/16.2MM</i>

HIP PLATE 120°

ABP-601	4	65MM
ABP-602	4	75MM
ABP-603	4	85MM

IN TITANIUM SAME SPECIFICATION IS FROM ABP-701 TO ABP-703





HIP PLATE 130°

ABP-801	4	65MM
ABP-802	4	70MM

IN TITANIUM SAME SPECIFICATIONS IS FROM ABP-901 TO ABP-902

ABP-1001	4	75MM
ABP-1002	4	80MM
ABP-1003	4	85MM
ABP-1004	4	90MM
ABP-1005	4	95MM
ABP-1006	4	100MM

IN TITANIUM SAME SPECIFICATIONS IS FROM ABP-2001 TO ABP-2006



PROFILE

<i>THICKNESS</i>	<i>5.7MM</i>
<i>WIDTH</i>	<i>16.2MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>
<i>U-PROFILE BLADE</i>	<i>6.6/16.2MM</i>

PROFILE

<i>THICKNESS</i>	<i>3.6MM</i>
<i>WIDTH</i>	<i>11.4MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>
<i>T-PROFILE BLADE</i>	<i>4.5/11.3MM</i>

CHILD HIP PLATES

CAT NO	HOLES	DISPLACEMENT	BLADE LENGTH
ABP-3001	3	8MM/80°	35MM
ABP-3002	3	8MM/90°	35MM
ABP-3003	3	8MM/100°	35MM

IN TITANIUM SAME SPECIFICATIONS IS FROM ABP-4001

TO ABP-4003





Ortho Innovations.

PAKISTAN

ADULT HIP PLATE 90°

ABP-5001 4 10MM 40MM
ABP-5002 4 10MM 60MM
IN TITANIUM SAME SPECIFICATIONS IS FORM ABP-6001
TO ABP-6002



PROFILE

<i>THICKNESS</i>	<i>5.7MM</i>
<i>WIDTH</i>	<i>16.12MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>
<i>U-PROFILE BLADE</i>	<i>6.5/1602MM</i>