




Ortho Innovations.

PAKISTAN

ORTHON BAISC INSTRUMENTS SET





Orthon Innovations.

PAKISTAN

QUICK COUPLING DRILL BIT

| | |
|---------|-------------|
| BIS-001 | 3.2MM/145MM |
| BIS-002 | 3.2MM/195MM |
| BIS-003 | 4.5MM/145MM |
| BSI-004 | 4.5MM/195MM |

JACOBS CHUCK DRILL BIT

| | |
|---------|-------------|
| BIS-005 | 3.2MM/130MM |
| BIS-006 | 3.2MM/180MM |
| BIS-007 | 4.5MM/130MM |


QUICK COUPLING TAP

| | |
|---------|-------------|
| BIS-008 | 4.5MM/125MM |
| BIS-009 | 4.5MM/175MM |
| BIS-010 | 6.5MM/195MM |

BIS-011
LARGE HEXOGONAL SCREW
DRIVER WITH HANDLE

BIS-012
HEXAGONAL SCREW
DRIVER WITH GROOVE

BIS-013
HOLDING SLEEVE
FOR BIS-005





Ortho Innovations.

PAKISTAN

BIS-014
COUNTER SINK



BIS-015
QUICK COUPLING TAP HANDLE



BIS-016
DOUBLE DRILL SLEEVE
4.5MM/3.2MM



BIS-017
DOUBLE DRILL SLEEVE
6.5MM/3.2MM



BIS-018
INSERT DRILL SLEEVE
4.5MM/3.2MM



BIS-019
DEPTH GAUGE





Ortho Innovations.

PAKISTAN

BIS-020
SHARP HOOK



BIS-021
COMBINATION WRENCH 11MM



BIS-022
NEUTRAL AND LOAD DRILL
SLEEVE 3.2MM



BIS-023
NEUTRAL AND LOAD
DRILL SLEEVE FOR LCP



BIS-024
UNIVERSAL DRILL SLEEVE
3.2MM/4.5MM



BIS-025
BENDING TEMPLATE 120MM
BIS-026
BENDING TEMPLATE 155MM
BIS-027
BENDING TEMPLATE 210MM





Ortho Innovations.

PAKISTAN

BIS-028
TENTION DEVICE 8MM



BIS-029
BENDING IRON



BIS-030
SOCKET WRENCH 11MM



BIS-031
WRENCH 8MM



BIS-032
QUICK COUPLING HEXAGONAL
SCREW DRIVER SHAFT





Ortho Innovations.

PAKISTAN

BIS-033
HEXAGONAL SCREW DRIVER
SHAFT WITH GROOVE



BIS-034
HEXAGONAL SCREW DRIVER
WITH UNIVERSAL JOINT

BIS-035
BENDING IRON PAIR



SELF CENTRING BONE HOLDING FORCEPS

BIS-036 SIZE 1:240MM LONG
BIS-037 SIZE2:260MM LONG
BIS-038 SIZE3:280MM LONG

BIS-039 LARGE REDUCTION FORCEPS
 POINTED 200MM LONG
BIS-040 REDUCTION FORCEPS
 WITH SERRATED JAWS 240MM LONG






Ortho Innovations.

PAKISTAN

ORTHON BASIC IMPLANTS SET



SEMI-TUBULAR PLATES

| CAT NO | HOLES | LENGTH |
|---------|-------|--------|
| DCP-001 | 2 | 39MM |
| DCP-002 | 3 | 55MM |
| DCP-003 | 4 | 71MM |
| DCP-004 | 5 | 87MM |
| DCP-005 | 6 | 103MM |
| DCP-006 | 7 | 119MM |
| DCP-007 | 8 | 135MM |
| DCP-008 | 9 | 151MM |
| DCP-009 | 10 | 167MM |
| DCP-010 | 11 | 183MM |
| DCP-011 | 12 | 199MM |

IN TITANIUM THE SAME SPECIFICATIONS ARE FROM DCP-101 TO DCP-111



| PROFILE | |
|--------------|-------------|
| THICKNESS | 1.0MM |
| WIDTH | 11.0MM |
| HOLE SPACING | 16.0&26.0MM |



NARROW DYNAMIC COMPRESSION PLATES

| | | |
|---------|----|-------|
| DCP-201 | 2 | 39MM |
| DCP-202 | 3 | 55MM |
| DCP-203 | 4 | 71MM |
| DCP-204 | 5 | 87MM |
| DCP-205 | 6 | 103MM |
| DCP-206 | 7 | 119MM |
| DCP-207 | 8 | 135MM |
| DCP-208 | 9 | 151MM |
| DCP-209 | 10 | 167MM |
| DCP-210 | 11 | 183MM |
| DCP-211 | 12 | 199MM |
| DCP-212 | 13 | 215MM |
| DCP-213 | 14 | 231MM |
| DCP-214 | 15 | 247MM |
| DCP-215 | 16 | 263MM |
| DCP-216 | 18 | 279MM |
| DCP-217 | 20 | 295MM |

IN TITANIUM SAME SPECIFICATIONS ARE IN DCP-301 TO DCP-317



| | |
|---------------------|--------|
| <i>PROFILE</i> | |
| <i>THICKNESS</i> | 4.7MM |
| <i>WIDTH</i> | 13.6MM |
| <i>HOLE SPACING</i> | 18.0MM |

BROAD DYNAMIC COMPRESSION PLATES

| | | |
|---------|----|-------|
| DCP-401 | 6 | 103MM |
| DCP-402 | 7 | 119MM |
| DCP-403 | 8 | 135MM |
| DCP-404 | 9 | 151MM |
| DCP-405 | 10 | 167MM |
| DCP-406 | 11 | 183MM |
| DCP-407 | 12 | 199MM |
| DCP-408 | 14 | 231MM |
| DCP-409 | 16 | 263MM |
| DCP-410 | 18 | 279MM |
| DCP-411 | 20 | 327MM |
| DCP-412 | 22 | 359MM |
| DCP-413 | 24 | 391MM |

IN TITANIUM SAME SPECIFICATIONS ARE IN DCP-501 TO DCP-513



| | |
|---------------------|---------------|
| <i>PROFILE</i> | |
| <i>THICKNESS</i> | 3.7MM |
| <i>WIDTH</i> | 12.0MM |
| <i>HOLE SPACING</i> | 16.0MM&26.0MM |



NARROW LIMITED CONTACT DYNAMIC COMPRESSION PLATES

| | | |
|----------|----|-------|
| NLCP-001 | 2 | 34MM |
| NLCP-002 | 3 | 52MM |
| NLCP-003 | 4 | 70MM |
| NLCP-004 | 5 | 88MM |
| NLCP-005 | 6 | 106MM |
| NLCP-006 | 7 | 124MM |
| NLCP-007 | 8 | 142MM |
| NLCP-008 | 9 | 160MM |
| NLCP-009 | 10 | 178MM |
| NLCP-010 | 11 | 196MM |
| NLCP-011 | 12 | 214MM |
| NLCP-012 | 14 | 250MM |
| NLCP-013 | 16 | 286MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM
NLCP-101 TO NLCP-113



PROFILE

| | |
|---------------------|---------------|
| <i>THICKNESS</i> | <i>4.7MM</i> |
| <i>WIDTH</i> | <i>13.6MM</i> |
| <i>HOLE SPACING</i> | <i>18.0MM</i> |

BROAD LIMITED CONTACT DYNAMIC COMPRESSION PLATES

| | | |
|----------|----|-------|
| BLCP-001 | 6 | 106MM |
| BLCP-002 | 7 | 124MM |
| BLCP-003 | 8 | 142MM |
| BLCP-004 | 9 | 160MM |
| BLCP-005 | 10 | 178MM |
| BLCP-006 | 11 | 196MM |
| BLCP-007 | 12 | 214MM |
| BLCP-008 | 13 | 232MM |
| BLCP-009 | 14 | 250MM |
| BLCP-010 | 15 | 268MM |
| BLCP-011 | 16 | 286MM |
| BLCP-012 | 17 | 304MM |
| BLCP-013 | 18 | 322MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM
BLCP-101 TO BLCP-113



PROFILE

| | |
|---------------------|---------------|
| <i>THICKNESS</i> | <i>6.1MM</i> |
| <i>WIDTH</i> | <i>17.6MM</i> |
| <i>HOLE SPACING</i> | <i>18.0MM</i> |



PAKISTAN

SPOON PLATES

| | | |
|---------|---|-------|
| DCP-601 | 5 | 100MM |
| DCP-602 | 6 | 120MM |
| DCP-603 | 7 | 140MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-701 TO DCP-703



PROFILE
THICKNESS 2.0MM
HEAD HOLE 33.0MM

T-PLATES

| | | |
|---------|----|-------|
| DCP-801 | 3 | 68MM |
| DCP-802 | 4 | 84MM |
| DCP-803 | 5 | 100MM |
| DCP-804 | 6 | 116MM |
| DCP-805 | 7 | 132MM |
| DCP-806 | 8 | 148MM |
| DCP-807 | 9 | 164MM |
| DCP-808 | 11 | 196MM |
| DCP-809 | 12 | 212MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-901 TO DCP-909



PROFILE
THICKNESS 2.0MM
WIDTH 17.0MM

T-BUTTRESS PLATES

| | | |
|----------|---|-------|
| DCP-1001 | 4 | 81MM |
| DCP-1002 | 5 | 96MM |
| DCP-1003 | 6 | 112MM |
| DCP-1004 | 7 | 128MM |
| DCP-1005 | 8 | 144MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-2001 TO DCP-2005



PROFILE
THICKNESS 2.0MM
WIDTH 17.0MM



L-BUTTRESS PLATES

LEFT ANGLED

| | | |
|----------|---|-------|
| DCP-3001 | 3 | 69MM |
| DCP-3002 | 4 | 89MM |
| DCP-3003 | 5 | 101MM |
| DCP-3004 | 6 | 117MM |


IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-4001 TO DCP-4004



RIGHT ANGLED

| | | |
|----------|---|-------|
| DCP-5001 | 3 | 69MM |
| DCP-5002 | 4 | 89MM |
| DCP-5003 | 5 | 101MM |
| DCP-5004 | 6 | 117MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-6001 TO DCP-6004



LATERAL TIBIAL HEAD BUTTRESS PLATES

LEFT ANGLED

| | | |
|----------|----|-------|
| DCP-7001 | 5 | 117MM |
| DCP-7002 | 7 | 149MM |
| DCP-7003 | 9 | 181MM |
| DCP-7004 | 11 | 213MM |
| DCP-7005 | 13 | 245MM |


IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-8001 TO DCP-8005



RIGHT ANGLED

| | | |
|----------|----|-------|
| DCP-9001 | 5 | 117MM |
| DCP-9002 | 7 | 149MM |
| DCP-9003 | 9 | 181MM |
| DCP-9004 | 11 | 213MM |
| DCP-9005 | 13 | 245MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM DCP-10001 TO DCP-10005





**CONDYLAR BUTTRESS PLATES
LEFT ANGLED**

| | | |
|-----------|----|-------|
| DCP-20001 | 7 | 157MM |
| DCP-20002 | 9 | 189MM |
| DCP-20003 | 11 | 221MM |
| DCP-20004 | 13 | 253MM |
| DCP-20005 | 15 | 285MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM
DCP-30001 TO DCP-30005



RIGHT ANGLED

| | | |
|-----------|----|-------|
| DCP-40001 | 7 | 157MM |
| DCP-40002 | 9 | 189MM |
| DCP-40003 | 11 | 221MM |
| DCP-40004 | 13 | 253MM |
| DCP-40005 | 15 | 285MM |

IN TITANIUM SAME SPECIFICATIONS ARE FROM
DCP-50001 TO DCP-50005



**DISTAL TIBIAL BUTTRESS PLATES
LEFT LEG**

| | | |
|-----------|----|-------|
| DCP-60001 | 7 | 122MM |
| DCP-60002 | 9 | 158MM |
| DCP-60003 | 11 | 194MM |

IN TITANIUM FROM DCP-70001 TO DCP-70003



RIGHT LEG

| | | |
|-----------|----|-------|
| DCP-80001 | 7 | 122MM |
| DCP-80002 | 9 | 158MM |
| DCP-80003 | 11 | 194MM |

IN TITANIUM FROM DCP-90001 TO DCP-90003





COBRA HEAD PLATES

| | | |
|------------|----|-------|
| DCP-100001 | 8 | 170MM |
| DCP-100002 | 9 | 186MM |
| DCP-100003 | 10 | 202MM |
| DCP-100004 | 11 | 218MM |

IN TITANIUM FROM DCP-200001 TO DCP-200004



| | |
|------------------|--------|
| <i>PROFILE</i> | |
| <i>THICKNESS</i> | 6.0MM |
| <i>WIDTH</i> | 16.0MM |

LCP PROXIMAL FEMORAL PLATE 4.5/5.0

FROM 2 HOLES TO 16 HOLES WITH DIFFERENCE OF TWO HOLES LEFT SIDE.

BLCF-001 TO BLCF-008

RIGHT SIDE

BLCF-001A TO BLCF-008A

LCP-PROXIMAL FEMUR HOOKS PLATE 4.5/5.0

FROM 2 HOLES TO 18 HOLES WITH DIFFERENCE OF TWO HOLES

BLCF-001B TO BLSP-009B





LCP BROAD 4.5/5.0 CURVED

| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| BLCP-001C | 12 | 229MM |
| BLCP-002C | 13 | 247MM |
| BLCP-003C | 14 | 265MM |
| BLCP-004C | 15 | 282MM |
| BLCP-005C | 16 | 300MM |
| BLCP-006C | 17 | 318MM |
| BLCP-007C | 18 | 336MM |
| BLCP-008C | 19 | 354MM |
| BLCP-009C | 20 | 372MM |
| BLCP-010C | 21 | 390MM |
| BLCP-011C | 22 | 408MM |
| BLCP-012C | 23 | 426MM |
| BLCP-013C | 24 | 443MM |
| BLCP-014C | 25 | 461MM |
| BLCP-015C | 26 | 479MM |

IN TITANIUM BLCP-101C TO BLCP-115C



LCP-PLT 4.5/5.0 PROXIMAL TIBIAL PLATE

| CAT NO | DIRECTION | HOLES | LENGTH |
|-----------|-----------|-------|--------|
| BLCP-001D | RIGHT | 5 | 140MM |
| BLCP-002D | LEFT | 5 | 140MM |
| BLCP-003D | RIFGHT | 7 | 180MM |
| BLCP-004D | LEFT | 7 | 180MM |
| BLCP-005D | RIGHT | 9 | 220MM |
| BLCP-006D | LEFT | 9 | 220MM |
| BLCP-007D | RIGHT | 11 | 260MM |
| BLCP-008D | LEFT | 11 | 260MM |
| BLCP-009D | RIGHT | 13 | 300MM |
| BLCP-010D | LEFT | 13 | 300MM |

IN TITANIUM BLCP-101D TO BLCP-110D





LCP MEDIAL PROXIMAL TIBIAL PLATE 3.5

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001E | RIGHT | 4 | 93MM |
| BLCP-002E | LEFT | 4 | 93MM |
| BLCP-003E | RIGHT | 6 | 119MM |
| BLCP-004E | LEFT | 6 | 119MM |
| BLCP-005E | RIGHT | 8 | 145MM |
| BLCP-006E | LEFT | 8 | 145MM |
| BLCP-007E | RIGHT | 10 | 171MM |
| BLCP-008E | LEFT | 10 | 171MM |
| BLCP-009E | RIGHT | 12 | 197MM |
| BLCP-010E | LEFT | 12 | 197MM |
| BLCP-011E | RIGHT | 14 | 223MM |
| BLCP-012E | LEFT | 14 | 223MM |
| BLCP-013E | RIGHT | 16 | 249MM |
| BLCP-014E | LEFT | 16 | 249MM |
| BLCP-015E | RIGHT | 18 | 275MM |
| BLCP-016E | LEFT | 18 | 275MM |
| BLCP-017E | RIGHT | 20 | 301MM |
| BLCP-018E | LEFT | 20 | 301MM |

IN TITANIUM BLCP-101E TO BLCP-118E



LCP MEDIAL PROXIMAL TIBIAL PLATE 4.5

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001F | RIGHT | 4 | 106MM |
| BLCP-002F | LEFT | 4 | 106MM |
| BLCP-003F | RIGHT | 6 | 142MM |
| BLCP-004F | LEFT | 6 | 142MM |
| BLCP-005F | RIGHT | 8 | 178MM |
| BLCP-006F | LEFT | 8 | 178MM |
| BLCP-007F | RIGHT | 10 | 214MM |
| BLCP-008F | LEFT | 10 | 214MM |
| BLCP-009F | RIGHT | 12 | 250MM |
| BLCP-010F | LEFT | 12 | 250MM |
| BLCP-011F | RIGHT | 14 | 286MM |
| BLCP-012F | LEFT | 14 | 286MM |
| BLCP-013F | RIGHT | 16 | 250MM |
| BLCP-014F | LEFT | 16 | 250MM |
| BLCP-015F | RIGHT | 18 | 286MM |





Ortho Innovations.

PAKISTAN

| | | | |
|-----------|-------|----|-------|
| BLCP-016F | LEFT | 18 | 286MM |
| BLCP-017F | RIGHT | 20 | 322MM |
| BLCP-018F | LEFT | 20 | 322MM |

IN TITANIUM BLCP-101F TO BLCP-118F



LCP-PROXIMAL TIBIAL PLATE 4.5/5.0 LATERAL

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001G | RIGHT | 4 | 82MM |
| BLCP-002G | LEFT | 4 | 82MM |
| BLCP-003G | RIGHT | 6 | 118MM |
| BLCP-004G | LEFT | 6 | 118MM |
| BLCP-005G | RIGHT | 8 | 154MM |
| BLCP-006G | LEFT | 8 | 154MM |
| BLCP-007G | RIGHT | 10 | 190MM |
| BLCP-008G | LEFT | 10 | 190MM |
| BLCP-009G | RIGHT | 12 | 226MM |
| BLCP-010G | LEFT | 12 | 226MM |
| BLCP-011G | RIGHT | 14 | 262MM |
| BLCP-012G | LEFT | 14 | 262MM |



LCP-DF-4.5/5.0

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001H | RIGHT | 5 | 156MM |
| BLCP-002H | LEFT | 5 | 156MM |
| BLCP-003H | RIGHT | 7 | 196MM |
| BLCP-004H | LEFT | 7 | 196MM |
| BLCP-005H | RIGHT | 9 | 236MM |
| BLCP-006H | LEFT | 9 | 236MM |
| BLCP-007H | RIGHT | 11 | 276MM |
| BLCP-008H | LEFT | 11 | 276MM |



Ortho Innovations.

PAKISTAN

BLCP-009H RIGHT 13 316MM
BLCP-010H LEFT 13 316MM

IN TITANIUM BLCP-101H TO BLCP-110H



LCP CONDYLAR PLATE 4.5/5.0

| CAT NO | DIRECTION | HOLES | LENGTGH | |
|-----------|-----------|-------|---------|--------|
| BLCP-001J | RIGHT | 6 | 170MM | |
| BLCP-002J | LEFT | 6 | 170MM | |
| BLCP-003J | RIGHT | 8 | 206MM | |
| BLCP-004J | LEFT | 8 | 206MM | |
| BLCP-005J | RIGHT | 10 | 242MM | CURVED |
| BLCP-006J | LEFT | 10 | 242MM | CURVED |
| BLCP-007J | RIGHT | 12 | 278MM | CURVED |
| BLCP-008J | LEFT | 12 | 278MM | CURVED |
| BLCP-009J | RIGHT | 14 | 314MM | CURVED |
| BLCP-010J | LEFT | 14 | 314MM | CURVED |
| BLCP-011J | RIGHT | 16 | 350MM | CURVED |
| BLCP-012J | LEFT | 16 | 350MM | CURVED |
| BLCP-013J | RIGHT | 18 | 386MM | CURVED |
| BLCP-014J | LEFT | 18 | 386MM | CURVED |
| BLCP-015J | RIGHT | 20 | 422MM | CURVED |
| BLCP-016J | LEFT | 20 | 422MM | CURVED |
| BLCP-017J | RIGHT | 22 | 458MM | CURVED |
| BLCP-018J | LEFT | 22 | 458MM | CURVED |



Ortho Innovations.

PAKISTAN



LCP 4.5/5.0 BROAD, CURVED

| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| BLCP-001K | 12 | 229MM |
| BLCP-002K | 13 | 247MM |
| BLCP-003K | 14 | 265MM |
| BLCP-004K | 15 | 282MM |
| BLCP-005K | 16 | 300MM |
| BLCP-006K | 17 | 318MM |
| BLCP-007K | 18 | 336MM |
| BLCP-008K | 19 | 354MM |
| BLCP-009K | 20 | 372MM |
| BLCP-010K | 21 | 390MM |
| BLCP-011K | 22 | 408MM |
| BLCP-012K | 23 | 426MM |
| BLCP-013K | 24 | 443MM |
| BLCP-014K | 25 | 461MM |
| BLCP-015K | 26 | 479MM |



Ortho Innovations.

PAKISTAN



LISS-PLT PLATE 5.0 TITANIUM

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001L | RIGHT | 5 | 141MM |
| BLCP-002L | LEFT | 5 | 141MM |
| BLCP-003L | RIGHT | 9 | 221MM |
| BLCP-004L | LEFT | 9 | 221MM |
| BLCP-005L | RIGHT | 13 | 301MM |
| BLCP-006L | LEFT | 13 | 301MM |



LISS-DF PLATE 5.0 TITANIUM

| CAT NO | DIRECTION | HOLES | LENGTGH |
|-----------|-----------|-------|---------|
| BLCP-001L | RIGHT | 5 | 156MM |
| BLCP-002L | LEFT | 5 | 156MM |
| BLCP-003L | RIGHT | 9 | 236MM |
| BLCP-004L | LEFT | 9 | 236MM |
| BLCP-005L | RIGHT | 13 | 316MM |
| BLCP-006L | LEFT | 13 | 316MM |

SCREW HOLE INSERT Ø 5.0MM

TITANIUM

HEXAGONAL SOCKET 3.5MMØ





LCP PLATE 4.5/5.0 NARROW

| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| NLCP-001A | 2 | 44MM |
| NLCP-002A | 3 | 62MM |
| NLCP-003A | 4 | 80MM |
| NLCP-004A | 5 | 98MM |
| NLCP-005A | 6 | 116MM |
| NLCP-006A | 7 | 134MM |
| NLCP-007A | 8 | 152MM |
| NLCP-008A | 9 | 170MM |
| NLCP-009A | 10 | 188MM |
| NLCP-010A | 11 | 206MM |
| NLCP-011A | 12 | 224MM |
| NLCP-012A | 13 | 242MM |
| NLCP-013A | 14 | 260MM |
| NLCP-014A | 15 | 278MM |
| NLCP-015A | 16 | 296MM |
| NLCP-016A | 18 | 332MM |
| NLCP-018A | 20 | 368MM |
| NLCP-019A | 22 | 404MM |
| NLCP-020A | 24 | 440MM |

IN TITANIUM NLCP-101A TO NLCP-120A



| PROFILE | |
|--------------|--------|
| THICKNESS | 4.6MM |
| WIDTH | 13.5MM |
| HOLE SPACING | 18.0MM |

LCP PLATE 4.5/5.0 BROAD

| | | |
|-----------|----|-------|
| NLCP-001B | 6 | 116MM |
| NLCP-002B | 7 | 134MM |
| NLCP-003B | 8 | 152MM |
| NLCP-004B | 9 | 170MM |
| NLCP-005B | 10 | 188MM |
| NLCP-006B | 11 | 206MM |
| NLCP-007B | 12 | 224MM |
| NLCP-008B | 13 | 242MM |
| NLCP-009B | 14 | 260MM |
| NLCP-010B | 15 | 278MM |
| NLCP-011B | 16 | 296MM |
| NLCP-012B | 18 | 332MM |
| NLCP-013B | 20 | 368MM |
| NLCP-014B | 22 | 404MM |
| NLCP-015B | 24 | 440MM |

IN TITANIUM NLCP-101B TO NLCP-115B



| PROFILE | |
|--------------|--------|
| THICKNESS | 5.2MM |
| WIDTH | 17.5MM |
| HOLE SPACING | 18.0MM |



LCP 4.0/5.0 BROAD CURVED



| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| NLCP-001C | 12 | 229MM |
| NLCP-002C | 13 | 247MM |
| NLCP-003C | 14 | 265MM |
| NLCP-004C | 15 | 2782M |
| NLCP-005C | 16 | 301MM |
| NLCP-006C | 18 | 337MM |

IN TITANIUM NLC-101C TO NLCP-106



LCP RECONSTRUCTION PLATE 4.5/5.0

| CAT NO | HOLES | LENGTH |
|----------|-------|--------|
| RLCP-001 | 3 | 56MM |
| RLCP-002 | 4 | 75MM |
| RLCP-003 | 5 | 94MM |
| RLCP-004 | 6 | 113MM |
| RLCP-005 | 7 | 132MM |
| RLCP-006 | 8 | 151MM |
| RLCP-007 | 9 | 170MM |
| RLCP-008 | 10 | 189MM |
| RLCP-009 | 11 | 208MM |
| RLCP-010 | 12 | 227MM |
| RLCP-011 | 13 | 246MM |
| RLCP-012 | 14 | 256MM |
| RLCP-013 | 15 | 284MM |
| RLCP-014 | 16 | 303MM |

IN TITANIUM RLCP-101 TO RLCP-114

PROFILE
THICKNESS 2.8MM STEEL
3.0MM TITANIUM
WIDTH 12.0MM
HOLE SPACING 19.0MM
DO NOT BEND OVER 15°



LCP-T-PLATE 4.5/5.0

| CAT NO | HOLES | LENGTH |
|----------|-------|--------|
| TLCP-001 | 3 | 68MM |
| TLCP-002 | 4 | 84MM |
| TLCP-003 | 5 | 100MM |
| TLCP-004 | 6 | 116MM |
| TLCP-005 | 7 | 134MM |
| TLCP-006 | 8 | 148MM |
| TLCP-007 | 10 | 180MM |
| TLCP-008 | 12 | 212MM |

IN TITANIUM TLCP-101 TO TLCP-108



LCP-L-BUTTRESS PLATE RIGHT, LEFT ANGLED

| CAT NO | HOLES | LENGTH |
|----------|-------|--------|
| LLCP-001 | 3 | 69MM |
| LLCP-002 | 4 | 85MM |
| LLCP-003 | 5 | 101MM |
| LLCP-004 | 6 | 117MM |

IN TITANIUM LLCP-101 TO LLCP-104



PROFILE
THICKNESS 2.0MM STEEL
2.5 MM TITANIUM
WIDTH 16.0MM

LCP -L-BUTTRESS PLATE 4.5/5.0 LEFT, RIGHT ANGLED

| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| LLCP-001A | 3 | 69MM |
| LLCP-002A | 4 | 85MM |
| LLCP-003A | 5 | 101MM |
| LLCP-004A | 6 | 117MM |

IN TITANIUM LLCP-101A TO LLCP-104A





Ortho Innovations.

PAKISTAN

LCP-T-BUTTRESS PLATE 4.5/5.0

| CAT NO | HOLES | LENGTH |
|-----------|-------|--------|
| TLCP-001A | 4 | 83MM |
| TLCP-002A | 5 | 98MM |
| TLCP-003A | 6 | 114MM |

IN TITANIUM TLCP=101A TO TLCP-103A



| PROFILE | |
|-----------|--------------------------------|
| THICKNESS | 2.0MM STEEL 2.5 MM TITANIUM |
| WIDTH | 17.0MM |

HOOK PLATE

HKP-001, 4, 5, 6 HOLES. LENGTH 123MM

