




Ortho Innovations.

PAKISTAN

**ORTHON CANNULATED SCREW
INSTRUMENTS 7.0MM**





Ortho Innovations.

PAKISTAN

CSL-001 SELF HOLDING SCREW FORCEPS



CSL-002 PARRELEL WIRE GUIDE FOR GUIDEWIRES 2.0MM



CSL-003 CANNULATED SCREW DRIVER



CSL-004 HOLDING SLEEVE FOR CSL-O-003




CSL-005 DIRECT MEASURING GUIDE



CSL-006 CANNULATED DRILL BIT 4.5MM/230MM



CSL-007 CANNULATED TAP 7.0MM



CSL-009 THREADED GUIDE WIRE 2.0MM/230MM





IMPLANTS (BOTH STAINLESS STEEL AND TITANIUM)

7.0MM CANNULATED SCREWS WITH 16MM THREADS

<i>CAT NO</i>	<i>LENGTH</i>
<i>CSL-101</i>	<i>30MM</i>
<i>CSL-102</i>	<i>35MM</i>
<i>CSL-103</i>	<i>40MM</i>
<i>CSL-104</i>	<i>45MM</i>
<i>CSL-105</i>	<i>50MM</i>
<i>CSL-106</i>	<i>55MM</i>
<i>CSL-107</i>	<i>60MM</i>
<i>CSL-108</i>	<i>65MM</i>
<i>CSL-109</i>	<i>70MM</i>
<i>CSL-110</i>	<i>75MM</i>
<i>CSL-111</i>	<i>80MM</i>
<i>CSL-112</i>	<i>85MM</i>
<i>CSL-113</i>	<i>90MM</i>
<i>CSL-114</i>	<i>95MM</i>
<i>CSL-115</i>	<i>100MM</i>
<i>CSL-116</i>	<i>105MM</i>
<i>CSL-117</i>	<i>110MM</i>
<i>CSL-118</i>	<i>115MM</i>
<i>CSL-119</i>	<i>120MM</i>



PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>SHAFT DIAMETER</i>	<i>4.5MM</i>
<i>THREAD DIAMETER</i>	<i>7.0MM</i>
<i>CORE DIAMETER</i>	<i>4.5MM</i>
<i>THREAD LENGTH</i>	<i>16.0MM</i>
<i>CANNULATION</i>	<i>2.1MM</i>
<i>DRILL BIT</i>	<i>4.5MM</i>
<i>TAP</i>	<i>7.0MM</i>

*IN TITANIUM SAME SPECIFICATIONS ARE FROM
CSL-201 TO CSL-219*



7.0MM CANNULATED SCREWS WITH 32MM THREADS

CSL-301	30MM
CSL-302	35MM
CSL-303	40MM
CSL-304	45MM
CSL-305	50MM
CSL-306	55MM
CSL-307	60MM
CSL-308	65MM
CSL-309	70MM
CSL-310	75MM
CSL-311	80MM
CSL-312	85MM
CSL-313	90MM
CSL-314	95MM
CSL-315	100MM
CSL-316	105MM
CSL-317	110MM
CSL-318	115MM
CSL-319	120MM



PROFILE	
HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
SHAFT DIAMETER	4.5MM
THREAD DIAMETER	7.0MM
CORE DIAMETER	4.5MM
THREAD LENGTH	32.0MM
CANNULATION	2.1MM
DRILL BIT	4.5MM
TAP	7.0MM

WASHER 13 MM DIAMETERS

CSL-320
IN TITANIUM SAME SPECIFICATIONS ARE
FROM CSL-401 TO CSL-420





Ortho Innovations.

PAKISTAN

ORTHON CANNULATED SCREW INSTRUMENTS 4.0MM



IN STERILIZING CONTAINER

ITEM NUMBER: INS-O-CSS

CAT NO	DESCRIPTION	QUANTITY
CSS-001	SELF HOLDING SCREW FORCEPS	1
CSS-002	DOUBLE DRILL SLEEVE 2.7/1.2MM	1
CSS-003	SMALL CANNULATED SCREW DRIVER	1
CSS-004	HOLDING SLEEVE FOR CSS-O-003	1
CSS-005	DIRECT MEASURING DEVICE	1
CSS-006	CANNULATED DRILL BIT 2.7MM	1
CSS-007	CANNULATED TAP	1
CSS-008	STYLET 2.0MM	1
CSS-009	THREADED GUIDE WIRE 2.0MM	1
CSS-010	DRILL SLEEVE 3.5/2.7MM WITH STOP	1
CSS-011	COUNTERSINK CANNULATED	1



Ortho Innovations.

PAKISTAN

CSS-001

QUICK COUPLING TAP HANDLE



CSS-002

DOUBLE DRILL SLEEVE 2.7/1.2MM



CSS-003

SELF HOLDING SCREW FORCEPS



CSS-004

CANULATED SCREW DRIVER



CSS-005

HOLDING SLEEVE FOR CSS-O-004



CSS-06

DIRECT MEASURING DEVICE





Ortho Innovations.

PAKISTAN

CSS-007

DRILL SLEEVE WITH STOP

3.5/2.7MM



CSS-008

CANNULATED DRILL BIT 2.7MM/160MM



CSS-009

CANNULATED TAP 4.0MM



CSS-010

CANNULATED COUNTERSINK



CSS-011

THREADED GUIDE PIN 1.25/150MM





Ortho Innovations.

PAKISTAN

IMPLANTS (IN BOTH STAINLESS STEEL AND TITANIUM)

4.0MM CANNULATED SCREWS SHORT THREADED

CAT NO	LENGTH
CSS-101	10MM
CSS-102	12MM
CSS-103	14MM
CSS-104	16MM
CSS-105	18MM
CSS-106	20MM
CSS-107	22MM
CSS-108	24MM
CSS-109	26MM
CSS-110	28MM
CSS-111	30MM
CSS-112	32MM
CSS-113	34MM
CSS-114	36MM
CSS-115	38MM
CSS-116	40MM
CSS-117	42MM
CSS-118	44MM
CSS-119	46MM
CSS-120	48MM
CSS-121	50MM



PROFILE

<i>HEAD DIAMETER</i>	<i>5.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>2.5MM</i>
<i>THREAD DIAMETER</i>	<i>4.0MM</i>
<i>CORE DIAMETER</i>	<i>2.5MM</i>
<i>SHAFT DIAMETER</i>	<i>2.7MM</i>
<i>DRILL BIT</i>	<i>2.7MM</i>
<i>TAP</i>	<i>4.0MM</i>

WASHER SMALL 7.0 MM DIAMETERS

CSS-122



IN TITANIUM SAME SPECIFICATIONS ARE FROM

CSS-201 TO CSS-221