




Ortho Innovations.

PAKISTAN

*ORTHON DYNAMIC HIP & DYNAMIC
CONDYLAR SCREW SET*





Ortho Innovations.

PAKISTAN

DHC-001 GUIDE WIRE 2.5MM/230MM

DHC-002 ANGLED GUIDE 135°



DHC-003 ANGLED GUIDE 150°



DHC-004 ANGLED GUIDE 95°



DHC-005 QUICK COUPLING T HANDLE



DHC-006 LAG SCREW WRENCH



DHC-007 CENTRINGSLEEVE FOR DHC-O-006





Ortho Innovations.

PAKISTAN

DHC-008

TAP 12.5MM ϕ



DHC-009

CENTRING SLEEVE FOR DHC-O-008



DHC-010

DIRECT MEASURING DEVICE



DHC-011

COUPLING SCREW FOR INSERTION



DHC-012

GUIDE SHAFT FOR COUPLING SCREW





Ortho Innovations.

PAKISTAN

DHC-013 COUPLING SCREW FOR REMOVE



DHC-014 PLATE IMPACTER



DHC-015 DHS TRIPPLE REAMER



DHC-016 DCS TRIPPLE REAMER



DHC-017 CANNULATED DRILL BIT 8.0MM





Ortho Innovations.

PAKISTAN

DHC-018

DHS REAMER



DHC-019

NUT FOR DHS REAMER



DHC-020

SHORT BARREL PLATE REAMER



DHC-021

DCS REAMER






Ortho Innovations.

PAKISTAN

ORTHON DHS IMPLANTS

BOTH IN STAINLESS STEEL AND TITANIUM



COMPRESSION SCREW

CAT NO	SIZE
DHS-001	36MM
IN TITANIUM DHS-101	



DYNAMIC HIP SCREWS

DHS-002	50MM
DHS-003	55MM
DHS-004	60MM
DHS-005	65MM
DHS-006	70MM
DHS-007	75MM
DHS-008	80MM
DHS-009	85MM
DHS-010	90MM
DHS-011	95MM
DHS-012	100MM
DHS-013	105MM
DHS-014	110MM
DHS-015	115MM
DHS-016	120MM
IN TITANIUM DHS-1002 TO DHS 116	






PAKISTAN

DYNAMIC HIP PLATES 135° LONG BARREL

CAT NO	HOLES	LENGTH
DHS-201	2	46MM
DHS-202	4	78MM
DHS-203	5	94MM
DHS-204	6	110MM
DHS-205	8	142MM
DHS-206	10	176MM
DHS-207	12	206MM
DHS-208	14	238MM
DHS-209	16	270MM

IN TITANIUM DHS-301 TO DHS-309



DYNAMIC HIP PLATES 135° SHORT BARREL

DHS-401	4	78MM
DHS-402	5	94MM
DHS-403	6	110MM

IN TITANIUM DHS-501 TO DHS-503



DYNAMIC HIP PLATES 150°

DHS-601	2	46MM
DHS-602	4	78MM
DHS-603	5	94MM
DHS-604	6	94MM
DHS-605	8	142MM
DHS-606	10	176MM
DHS-607	12	206MM



Ortho Innovations.

PAKISTAN

DHS-608 14 238MM
DHS-609 16 270MM
IN TITANIUM DHS-701 TO DHS-709


DYNAMIC HIP PLATES 130°

DHS-801 4 78MM
DHS-802 5 94MM
DHS-803 6 110MM
IN TITANIUM DHS-901 TO DHS- 903



DYNAMIC CONDYLAR PLATES 95°

DCS-001 6 100MM
DCS-002 8 130MM
DCS-003 10 163MM
DCS-004 12 198MM
DCS-005 14 225MM
DCS-006 16 260MM
IN TITANIUM DCS-101 TO DCS-106





LCP DHS PLATE 130° LONG BARREL (38MM)

CAT NO	HOLES	LENGTH
LDHS-001	2	60MM
LDHS-002	3	76MM
LDHS-003	4	92MM
LDHS-004	5	108MM
LDHS-005	6	124MM
LDHS-006	8	156MM
LDHS-007	10	188MM
LDHS-008	12	220MM
LDHS-009	14	252MM

IN TITANIUM LDHS001A TO LDHS-009A

LCP DHS PLATE 135 ° LONG BARREL (38MM)

CAT NO	HOLES	LENGTH
LDHS-201	2	60MM
LDHS-202	3	76MM
LDHS-203	4	92MM
LDHS-204	5	108MM
LDHS-205	6	124MM
LDHS-206	8	156MM
LDHS-207	10	188MM
LDHS-208	12	220MM
LDHS-209	14	252MM
LDHS-210	16	284MM
LDHS-211	18	316MM
LDHS-212	20	348MM

IN TITANIUM LDHS201A TO LDHS-212A

LCP DHS PLATE 140 ° LONG BARREL (38MM)

CAT NO	HOLES	LENGTH
LDHS-301	2	60MM
LDHS-302	3	76MM
LDHS-303	4	92MM
LDHS-304	5	108MM
LDHS-305	6	124MM
LDHS-306	8	156MM
LDHS-307	10	188MM
LDHS-308	12	220MM
LDHS-309	14	252MM





Ortho Innovations.

PAKISTAN

LDHS-310 16 284MM
IN TITANIUM LDHS301A TO LDHS-310A
LCP DHS PLATE 145 ° LONG BARREL (38MM)

CAT NO	HOLES	LENGTH
LDHS-401	2	60MM
LDHS-402	3	76MM
LDHS-403	4	92MM
LDHS-404	5	108MM
LDHS-405	6	124MM
LDHS-406	8	156MM
LDHS-407	10	188MM
LDHS-408	12	220MM
LDHS-409	14	252MM
LDHS-410	16	284MM
LDHS-411	18	316MM
LDHS-412	20	348MM

IN TITANIUM LDHS401A TO LDHS-412A
LCP DHS PLATE 150 ° LONG BARREL (38MM)

CAT NO	HOLES	LENGTH
LDHS-501	2	60MM
LDHS-502	3	76MM
LDHS-503	4	92MM
LDHS-504	5	108MM
LDHS-505	6	124MM
LDHS-506	8	156MM
LDHS-507	10	188MM
LDHS-508	12	220MM
LDHS-509	14	252MM
LDHS-510	16	284MM
LDHS-511	18	316MM
LDHS-512	20	348MM

IN TITANIUM LDHS501A TO LDHS-512A



Ortho Innovations.

PAKISTAN

LCP DHS PLATE 130 ° SHORT BARREL (25MM)

CAT NO	HOLES	LENGTH
LDHS-601	2	60MM
LDHS-602	4	92MM

IN TITANIUM LDHS-601A TO LDHS-602A

LCP DHS PLATE 135 ° SHORT BARREL (25MM)

CAT NO	HOLES	LENGTH
LDHS-701	2	60MM
LDHS-702	4	92MM
LDHS-703	5	108MM
LDHS-704	6	124MM

IN TITANIUM LDHS-701A TO LDHS-704A

LCP DHS PLATE 140 ° SHORT BARREL (25MM)

CAT NO	HOLES	LENGTH
LDHS-801	2	60MM
LDHS-802	4	92MM

IN TITANIUM LDHS-801A TO LDHS-802A

LCP DHS PLATE 145 ° SHORT BARREL (25MM)

CAT NO	HOLES	LENGTH
LDHS-901	2	60MM
LDHS-902	4	92MM

IN TITANIUM LDHS-901A TO LDHS-902A

LCP DHS PLATE 150 ° SHORT BARREL (25MM)

CAT NO	HOLES	LENGTH
LDHS-1001	2	60MM
LDHS-1002	4	92MM

IN TITANIUM LDHS-1001A TO LDHS-1002A





**LCP DHS PLATE WITH COLLER 135° LONG BARREL (38MM)
TITANIUM**

CAT NO	HOLES	LENGTH
LDHSC-001	3	96MM
LDHSC-002	4	112MM
LDHSC-003	5	128MM



**LCP DHS PLATE WITH COLLER 140° LONG BARREL (38MM)
TITANIUM**

CAT NO	HOLES	LENGTH
LDHSC-001A	3	96MM
LDHSC-002A	4	112MM
LDHSC-003A	5	128MM

DHS TROCHANTER STABILIZING PLATE 138 MM LENGTHS

STEEL	TITANIUM
DHST-001	DHST-001A



DHS TROCHANTER STABILIZING PLATE 148 MM LENGTHS

STEEL	TITANIUM
DHST-002	DHST-002A





Ortho Innovations.

PAKISTAN

DHS TROCHANTER STABILIZING PLATE 130 MM LENGTHS LOCKING

STEEL

TITANIUM

DHST-003

DHST-003A

