




Ortho Innovations.

PAKISTAN

ORTHON GENERAL INSTRUMENTS SET



IN STERILIZING CONTAINER

ITEM NUMBER: INS-O-GIS

CAT NO	DESCRIPTION	QUANTITY
RETRACTERS		
GIS-001	RETRACTER 8MM/220MM SHORT NARROW TIP	2
GIS-002	RETRACTER 18MM/235MM SHORT NARROW TIP	2
GIS-003	RETRACTER 43MM/240MM NARROW TIP	1
GIS-004	RETRACTER 70MM/240MM NARROW TIP	1
GIS-005	RETRACTER 43MM/235MM LONG ANGLED BLADE	1
GIS-006	RETRACTER 22MM/250MM WIDE TIP	1
GIS-007	RETRACTER 35MM/300MM WIDE TIP	1
GIS-008	RETRACTER 18MM/235MM LONG NARROW TIP	1
GIS-009	RETRACTER 24MM/270MM LONG WIDE TIP	1
GIS-010	SMALL RETRACTER 6MM/160MM SHORT NARROW TIP	1
GIS-011	SMALL RETRACTER 8MM/160MM SHORT NARROW TIP	1
GIS-012	SMALL RETRACTER 15MM	1
PERIOSTEAL ELEVATOR		
GIS-013	PERIOSTEAL ELEVATOR 3/200MM, CURVED BLADE, STRAIGHT EDGE	1
GIS-014	PERIOSTEAL ELEVATOR 6/200MM CURVED BLADE, STRAIGHT EDGE	



		PAKISTAN
GIS-015	PERIOSTEAL ELEVATOR 6/200MM, CURVED BALDE, ROUND EDGE	1 1
GIS-016	PERIOSTEAL ELEVATOR 13MM/200MM STRAIGHT BLADE STRAIGHT EDGE	1
GIS-017	PERIOSTEAL ELEVATOR 13MM/200MM CURVED BLADE STRAIGHT EDGE	1
GIS-018	PERIOSTEAL ELEVATOR 14/200MM STRAIGHT BLADE ROUND EDGE	1
GIS-019	PERIOSTEAL ELEVATOR 14/200MM CURVED BLADE ROUND EDGE	1
GIS-020	PERIOSTEAL ELEVATOR 20/200MM STRAIGHT BLADE ROUND EDGE	1
GOUGES		
GIS-021	GOUGE 5MM STRAIGHT	1
GIS-022	GOUGE 10MM STRAIGHT	1
GIS-023	GOUGE 15MM STRAIGHT	1
GIS-024	GOUGE 5MM CURVED	1
GIS-025	GOUGE 10MM CURVED	1
GIS-026	GOUGE 15MM CURVED	1
CHIESELS		
GIS-027	CHIESEL HANDLE	1
GIS-028	CHIESEL BLADE 5MM	1
GIS-029	CHIESEL BLADE 10MM	1
GIS-030	CHIESEL BLADE 16MM	1
GIS-031	CHIESEL BLADE 25MM	1
OSTEOTOMES		
GIS-032	4MM/STRAIGHT	1
GIS-033	6MM/STRAIGHT	1
GIS-034	8MM/STRAIGHT	1
GIS-035	10MM/STRAIGHT	1
GIS-036	15MM/STRAIGHT	1
GIS-037	20MM/STRAIGHT	1
GIS-038	25MM/STRAIGHT	1
GIS-039	30MM/STRAIGHT	1



GOUGES

GIS-040	4MM/STRAIGHT	1
GIS-041	6MM/STRAIGHT	1
GIS-042	8MM/STRAIGHT	1
GIS-043	10MM/STRAIGHT	1
GIS-044	15MM/STRAIGHT	1
GIS-045	20MM/STRAIGHT	1
GIS-046	25MM/STRAIGHT	1
GIS-047	30MM/STRAIGHT	1

CHISELS

GIS-048	4MM/STRAIGHT	1
GIS-049	6MM/STRAIGHT	1
GIS-050	8MM/STRAIGHT	1
GIS-051	10MM/STRAIGHT	1
GIS-052	15MM/STRAIGHT	1
GIS-053	20MM/STRAIGHT	1
GIS-054	25MM/STRAIGHT	1
GIS-055	30MM/STRAIGHT	1
GIS-056	10MM/CURVED	1
GIS-057	16MM/CURVED	1

BONE HOOKS

GIS-057	SMALL SHARP	1
GIS-058	SMALL BLUNT	1
GIS-059	LARGE SHARP	1
GIS-060	LARGE BLUNT	1

BONE HOOKS WITH T-HANDLE

GIS-061	LARGE SHARP	1
---------	-------------	---



Orthon Innovations.

PAKISTAN

HAMMERS

GIS-062	150G	1
GIS-063	350G	1
GIS-064	500G	1
GIS-065	700G	1
GIS-066	900G	1
GIS-067	NYLON FACED HAMMER	1
GIS-068	FIBRE MALLET	1



Ortho Innovations.

PAKISTAN

GIS-001 SHORT NARROW TIP RETRACTER 8MM/220MM



GIS-002 SHORT NARROW TIP RETRACTER 18MM/235MM



GIS-003 NARROW TIP RETRACTER 43MM/240MM



GIS-004 NARROW TIP RETRACTER 70MM/240MM



GIS-005 LONG ANGLED BLADE RETRACTER 43MM/235MM



GIS-006 WIDE TIP RETRACTER 22MM/250MM





Ortho Innovations.

PAKISTAN

GIS-007

WIDE TIP RETRACTER 35M/300MM



GIS-008

LONG NARROW TIP RETARCTER 18MM/235M



GIS-009

LONG WIDE TIP RETRACTER 24MM/270MM



GIS-010

SMALL SHORT NARROW TIP R ETRACTER 6MM/160MM



GIS-011

SMALL SHORT NARROW TIP RETRACTER 8MM/160MM



GIS-012

SMALL RETRACTER 15MM





Ortho Innovations.

PAKISTAN

GIS-013

PERIOSTEAL ELEVATOR, CURVED BLADE, STRAIGHT EDGE 3MM/ 200MM



GIS-014

PERIOSTEAL ELEVATOR, CURVED BLADE, STRAIGHT EDGE 6MM/200MM



GIS-015

PERIOSTEAL ELEVATOR, CURVED BLADE, ROUND EDGE 6MM/200MM



GIS-016

PERIOSTEAL ELEVATOR, STRAIGHT BLADE, STRAIGHT EDGE, 13MM/200MM



GIS-017

PERIOSTEAL ELEVATOR, CURVED BLADE, STRAIGHT EDGE, 13MM/200MM





Ortho Innovations.

PAKISTAN

GIS-018

PERIOSTEAL ELEVATOR, STRAIGHT BLADE, ROUND EDGE, 14MM/200MM



GIS-019

PERIOSTEAL ELEVATOR, CURVED BLADE, ROUND EDGE, 14MM/200MM



GIS-020

PERIOSTEAL ELEVATOR, STRAIGHT BLADE, ROUND EDGE, 20MM/200MM



GOUGES

GIS-021

5MM WIDE STRAIGHT

GIS-022

10MM WIDE STRAIGHT

GIS-023

15MM WIDE STRAIGHT

GIS-024

5MM WIDE CURVED

GIS-025

10MM WIDE CURVED

GIS-026

15MM WIDE CURVED





Ortho Innovations.

PAKISTAN

CHIESELS

GIS-027
GIS-028
GIS-029
GIS-030
GIS-031

CHIESEL HANDLE
CHIESEL BLADE 5MM
CHIESEL BLADE 10MM
CHIESEL BLADE 15MM
CHIESEL BLADE 25MM



OSTEOTOMES

GIS-032
GIS-033
GIS-034
GIS-035
GIS-036
GIS-037
GIS-038
GIS-039

4MM WIDE STRAIGHT
6 MM WIDE STRAIGHT
8 MM WIDE STRAIGHT
10 MM WIDE STRAIGHT
15 MM WIDE STRAIGHT
20 MM WIDE STRAIGHT
25 MM WIDE STRAIGHT
30 MM WIDE STRAIGHT



GOUGES

GIS-040
GIS-041
GIS-042
GIS-043
GIS-044
GIS-045
GIS-046
GIS-047

4MM WIDE STRAIGHT
6 MM WIDE STRAIGHT
8 MM WIDE STRAIGHT
10 MM WIDE STRAIGHT
15 MM WIDE STRAIGHT
20 MM WIDE STRAIGHT
25 MM WIDE STRAIGHT
30 MM WIDE STRAIGHT





CHIESELS

GIS-048	4MM WIDE STRAIGHT
GIS-049	6 MM WIDE STRAIGHT
GIS-050	8 MM WIDE STRAIGHT
GIS-051	10 MM WIDE STRAIGHT
GIS-052	15 MM WIDE STRAIGHT
GIS-053	20 MM WIDE STRAIGHT
GIS-054	25 MM WIDE STRAIGHT
GIS-055	30 MM WIDE STRAIGHT
GIS-056	10 MM WIDE CURVED
GIS-057	16 MM WIDE CURVED



BONE HOOKS

GIS-058	SMALL SHARP
GIS-059	SMALL BLUNT
GIS-060	LARGE SHARP
GIS-061	LARGE BLUNT



GIS-062	LARGE SHARP WITH T HANDLE
---------	---------------------------





Orthon Innovations.

PAKISTAN

HAMMERS

GIS-063 WEIGHT 150G
GIS-064 WEIGHT 350G
GIS-065 WEIGHT 500G
GIS-066 WEIGHT 700G
GIS-067 WEIGHT 900G



GIS-068 NYLON FACED



GIS-069 FIBRE MALLET






Ortho Innovations.

PAKISTAN

ORTHON BONE FORCEPS SET



IN STERILIZING CONTAINER

ITEM NUMBER: INS-O-BFS

CAT NO	DESCRIPTION	QUANTITY
BFS-001	SELF CENTRIN BONE HOLDING SIZE 0 LENGTH 190MM	1
BFS-002	SELF CENTRING BONE FORCEPS SIZE 1 LENGTH 240MM	1
BFS-003	SELF CENTRING BONE FORCEPS SIZE 2 LENGTH 260MM	1
BFS-004	SELF CENTRING BONE FORCEPS SIZE 3 LENGTH 280MM	2
BFS-005	SMALL REDUCTION FORCEPS WITH POINTS-RATCHET LOCK 130MM LENGTH	1
BFS-006	LARGE REDUCTION FORCEPS, WITH POINTS-RATCHET LOCK 200MM LENGTH	2
BFS-007	REDUCTION FORCEPS, SERRATED JAWS EXTRA SMALL RATCHETS LOCK 140MM LENGTH	
BFS-008	REDUCTION FORCEPS SERRATED JAWS, EXTRA SMALL THREADED LOCK 170MM LENGTH	1
BFS-009	LARGE REDUCTION FORCEPS SERRATED JAWS, THREADED LOCK 240MMLENGTH	1
BFS-010	REDUCTION FORCEPS WITH POINTS-RATCHET LOCK 165MM LENGTH	1



Ortho Innovations.

PAKISTAN

BFS-001

SELF CENTRING BONE FORCEPS 190MM LENGTH




BFS-002

SELF CENTRING BONE FORCEPS 240MM LENGTH



BFS-003

SELF CENTRING BONE HOLDING FORCEPS 260MM LENGTH



BFS-004

SELF CENTRING BONE FORCEPS 280MM LENGTH





PAKISTAN

BFS-005

SMALL REDUCTION FORCEPS, WITH POINTS,
RACHET LOCK 130MM LENGTH



BFS-006

LARGE REDUCTION FORCEPS WITH POINTS
RACHET LOCK, 200MM LENGTH



BFS-007

REDUCTION FORCEPS SERRATED JAWS,
EXTRA SMALL RACHET LOCK, 140MM LENGTH



BFS-008

SMALL REDUCTION FORCEPS,
SERRATED JAWS, THREADED LOCK, 170MM LENGTH



BFS-009

LARGE REDUCTION FORCEPS, SERRATED JAWS,
THREADED LOCK, 240MM LENGTH



BFS-010

REDUCTION FORCEPS, WITH POINTS
RATCHET LOCK, 165MM LENGTH

