




Ortho Innovations.

PAKISTAN



ORTHON LARGE SCREW SET

4.5MM CORTEX SCREWS

LSS-001	14MM
LSS-002	16MM
LSS-003	18MM
LSS-004	20MM
LSS-005	22MM
LSS-006	24MM
LSS-007	26MM
LSS-008	28MM
LSS-009	30MM
LSS-010	32MM
LSS-011	34MM
LSS-012	36MM
LSS-013	38MM

PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>THREAD DIAMETER</i>	<i>4.5MM</i>
<i>CORE DIAMETER</i>	<i>3.0MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>4.5MM</i>



Ortho Innovations.

PAKISTAN

LSS-014	40MM
LSS-015	42MM
LSS-016	44MM
LSS-017	46MM
LSS-018	48MM
LSS-019	50MM
LSS-020	52MM
LSS-021	54MM
LSS-022	56MM
LSS-023	58MM
LSS-024	60MM
LSS-025	62MM
LSS-026	64MM
LSS-027	66MM
LSS-028	68MM
LSS-029	70MM
LSS-030	72MM
LSS-031	74MM
LSS-032	76MM
LSS-033	78MM
LSS-034	80MM
LSS-035	82MM
LSS-036	84MM
LSS-037	86MM
LSS-038	88MM
LSS-039	90MM
LSS-040	95MM
LSS-041	100MM
LSS-042	105MM
LSS-043	110MM
LSS-044	115MM
LSS-045	120MM

**IN TITANIUM LSS-201 TO LSS-245
4.5MM MALLEOLAR SCREW**

LSS-101	25MM
LSS-102	30MM
LSS-103	35MM
LSS-104	40MM
LSS-105	45MM



PROFILE

HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	4.5MM
CORE DIAMETER	3.0MM



Ortho Innovations.

PAKISTAN

LSS-106	50MM
LSS-107	55MM
LSS-108	60MM
LSS-109	65MM
LSS-110	70MM

IN TITANIUM FROM LSS-301 TO LSS-310

4.5MM CORTEX SCREWS SELF TAPPING

LSS-401	14MM
LSS-402	16MM
LSS-403	18MM
LSS-404	20MM
LSS-405	22MM
LSS-406	24MM
LSS-407	26MM
LSS-408	28MM
LSS-409	30MM
LSS-410	32MM
LSS-411	34MM
LSS-412	36MM
LSS-413	38MM
LSS-414	40MM
LSS-415	42MM
LSS-416	44MM
LSS-417	46MM
LSS-418	48MM
LSS-419	50MM
LSS-420	52MM
LSS-421	54MM
LSS-422	56MM
LSS-423	58MM
LSS-424	60MM
LSS-425	62MM
LSS-426	64MM
LSS-427	66MM
LSS-428	68MM
LSS-429	70MM





Ortho Innovations.

PAKISTAN

LSS-430	72MM
LSS-431	74MM
LSS-432	76MM
LSS-433	78MM
LSS-434	80MM
LSS-435	82MM
LSS-436	84MM
LSS-437	86MM
LSS-438	88MM
LSS-439	90MM
LSS-440	95MM
LSS-441	100MM
LSS-442	105MM
LSS-443	110MM
LSS-444	115MM
LSS-445	120MM

IN TITANIUM LSS-501 TO LSS-545

6.5MM CANCELLOUS SCREW 16MM THREADS

LSS-601	25MM
LSS-602	30MM
LSS-603	35MM
LSS-604	40MM
LSS-605	45MM
LSS-606	50MM
LSS-607	55MM
LSS-608	60MM
LSS-609	65MM
LSS-610	70MM
LSS-611	75MM
LSS-612	80MM
LSS-613	85MM
LSS-614	90MM
LSS-615	95MM
LSS-616	100MM
LSS-617	105MM
LSS-618	110MM
LSS-619	115MM
LSS-620	120MM

IN TITANIUM FROM LSS-701 TO LSS-720

PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>THREAD DIAMETER</i>	<i>4.5MM</i>
<i>CORE DIAMETER</i>	<i>3.0MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>



PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>SHAFT DIAMETER</i>	<i>4.5MM</i>
<i>THREAD DIAMETER</i>	<i>6.5MM</i>
<i>CORE DIAMETER</i>	<i>3.0MM</i>
<i>THREAD LENGTH</i>	<i>16MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>6.5MM</i>



6.5MM CANCELLOUS SCREW 32MM THREADS

LSS-801	25MM
LSS-802	30MM
LSS-803	35MM
LSS-804	40MM
LSS-805	45MM
LSS-806	50MM
LSS-807	55MM
LSS-808	60MM
LSS-809	65MM
LSS-810	70MM
LSS-811	75MM
LSS-812	80MM
LSS-813	85MM
LSS-814	90MM
LSS-815	95MM
LSS-816	100MM
LSS-817	105MM
LSS-818	110MM
LSS-819	115MM
LSS-820	120MM



PROFILE	
HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
SHAFT DIAMETER	4.5MM
THREAD DIAMETER	6.5MM
CORE DIAMETER	3.0MM
THREAD LENGTH	32MM
DRILLBIT	3.2MM
TAP	6.5MM

IN TITANIUM LSS-901 TO LSS-920

6.5MM CANCELLOUS SCREW FULL THREAD

LSS-1001	25MM
LSS-1002	30MM
LSS-1003	35MM
LSS-1004	40MM
LSS-1005	45MM
LSS-1006	50MM
LSS-1007	55MM
LSS-1008	60MM
LSS-1009	65MM
LSS-1010	70MM
LSS-1011	75MM
LSS-1012	80MM
LSS-1013	85MM
LSS-1014	90MM
LSS-1015	95MM



PROFILE	
HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	6.5MM
DRILL BIT	3.2MM
TAP	6.5MM



Ortho Innovations.

PAKISTAN

LSS-1016 100MM
LSS-1017 105MM
LSS-1018 110MM
LSS-1019 115MM
LSS-1020 120MM
IN TITANIUM LSS-1101 TO LSS-1120

LARGE WASHER 13 MM DIAMETERS

LSS-1201



SCREW HOLDING FORCEPS
LSS-1202





Ortho Innovations.

PAKISTAN

LOCKING SCREWS

5.0MM Ø SELF TAPPING

LSS-001A	14MM
LSS-002A	16MM
LSS-003A	18MM
LSS-004A	20MM
LSS-005A	22MM
LSS-006A	24MM
LSS-007A	26MM
LSS-008A	28MM
LSS-009A	30MM
LSS-010A	32MM
LSS-011A	34MM
LSS-012A	36MM
LSS-013A	38MM
LSS-014A	40MM
LSS-015A	42MM
LSS-016A	44MM
LSS-017A	46MM
LSS-018A	48MM
LSS-019A	50MM
LSS-020A	52MM
LSS-021A	54MM
LSS-022A	56MM
LSS-023A	58MM
LSS-024A	60MM
LSS-025A	62MM
LSS-026A	64MM
LSS-027A	66MM
LSS-028A	68MM
LSS-029A	70MM
LSS-030A	72MM
LSS-031A	74MM
LSS-032A	76MM
LSS-033A	78MM
LSS-034A	80MM
LSS-035A	82MM
LSS-036A	84MM

PROFILE

HEAD DIAMETER	6.6MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	5.0MM
CORE DIAMETER	4.4MM
DRILL BIT	3.2MM
TAP	4.5MM





Ortho Innovations.

PAKISTAN

LSS-037 A	86MM
LSS-038 A	88MM
LSS-039 A	90MM

IN TITANIUM LSS-101A TO LSS-139A

LOCKING SCREW STARDRIVE Ø 5.0 SELF TAPPING

LSS-001B	14MM
LSS-002 B	16MM
LSS-003 B	18MM
LSS-004 B	20MM
LSS-005 B	22MM
LSS-006 B	24MM
LSS-007 B	26MM
LSS-008 B	28MM
LSS-009 B	30MM
LSS-010 B	32MM
LSS-011 B	34MM
LSS-012 B	36MM
LSS-013 B	38MM
LSS-014 B	40MM
LSS-015 B	42MM
LSS-016 B	44MM
LSS-017 B	46MM
LSS-018 B	48MM
LSS-019 B	50MM
LSS-020 B	52MM
LSS-021 B	54MM
LSS-022 B	56MM
LSS-023 B	58MM
LSS-024 B	60MM
LSS-025 B	62MM
LSS-026 B	64MM
LSS-027 B	66MM
LSS-028 B	68MM
LSS-029 B	70MM
LSS-030 B	72MM
LSS-031 B	74MM
LSS-032 B	76MM
LSS-033 B	78MM
LSS-034 B	80MM

PROFILE

<i>HEAD DIAMETER</i>	<i>6.6MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>THREAD DIAMETER</i>	<i>5.0MM</i>
<i>CORE DIAMETER</i>	<i>4.4MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>4.5MM</i>





Ortho Innovations.

PAKISTAN

LSS-035 B	82MM
LSS-036 B	84MM
LSS-037 B	86MM
LSS-038 B	88MM
LSS-039 B	90MM

IN TITANIUM LSS-101B TO LSS-139B

LOCKING SCREW Ø 5.0 WELF DRILLING

LSS-001B	14MM
LSS-002 B	16MM
LSS-003 B	18MM
LSS-004 B	20MM
LSS-005 B	22MM
LSS-006 B	24MM
LSS-007 B	26MM
LSS-008 B	28MM
LSS-009 B	30MM
LSS-010 B	32MM
LSS-011 B	34MM
LSS-012 B	36MM
LSS-013 B	38MM
LSS-014 B	40MM
LSS-015 B	42MM
LSS-016 B	44MM
LSS-017 B	46MM
LSS-018 B	48MM
LSS-019 B	50MM
LSS-020 B	52MM
LSS-021 B	54MM
LSS-022 B	56MM
LSS-023 B	58MM
LSS-024 B	60MM
LSS-025 B	62MM
LSS-026 B	64MM
LSS-027 B	66MM
LSS-028 B	68MM
LSS-029 B	70MM
LSS-030 B	72MM
LSS-031 B	74MM

PROFILE

<i>HEAD DIAMETER</i>	<i>6.6MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>THREAD DIAMETER</i>	<i>5.0MM</i>
<i>CORE DIAMETER</i>	<i>4.4MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>4.5MM</i>





LSS-032 B	76MM
LSS-033 B	78MM
LSS-034 B	80MM

IN TITANIUM LSS-101B TO LSS-134B

LOCKING SCREW STARDRIVE Ø 5.0 SELF DRILLING

LSS-001C	14MM
LSS-002C	16MM
LSS-003C	18MM
LSS-004C	20MM
LSS-005C	22MM
LSS-006C	24MM
LSS-007C	26MM
LSS-008C	28MM
LSS-009C	30MM
LSS-010C	32MM
LSS-011C	34MM
LSS-012C	36MM
LSS-013C	38MM
LSS-014C	40MM
LSS-015C	42MM
LSS-016C	44MM
LSS-017C	46MM
LSS-018C	48MM
LSS-019C	50MM
LSS-020C	52MM
LSS-021C	54MM
LSS-022C	56MM
LSS-023C	58MM
LSS-024C	60MM
LSS-025C	62MM
LSS-026C	64MM
LSS-027C	66MM
LSS-028C	68MM
LSS-029C	70MM
LSS-030C	72MM
LSS-031C	74MM
LSS-032C	76MM
LSS-033C	78MM
LSS-034C	80MM

IN TITANIUM LSS-101C TO LSS-134C

PROFILE	
HEAD DIAMETER	6.6MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	5.0MM
CORE DIAMETER	4.4MM
DRILL BIT	3.2MM
TAP	4.5MM





Ortho Innovations.

PAKISTAN

LOCKING SCREW Ø 5.0 MM SELF TAPPING FOR PERIPROSTHETIC FRACTURES

CAT NO	LENGTH
LSS-001P	8MM
LSS-002P	10MM
LSS-003P	12MM
LSS-004P	14MM
LSS-005P	16MM
LSS-006P	18MM

IN TITANIUM LSS101P TO LSS-106P



LOCKING SCREW Ø 5.0 MM SELF TAPPING FOR PERIPROSTHETIC FRACTURES ATADRIVE

CAT NO	LENGTH
LSS-001PS	8MM
LSS-002PS	10MM
LSS-003PS	12MM
LSS-004PS	14MM
LSS-005PS	16MM
LSS-006PS	18MM

IN TITANIUM LSS-101PS TO LSS-106PS



LCP SPACER Ø 5.0MM SIZE 2MM

LSS-001SP	STEEL
LSS-101SP	TITANIUM



**SCREW HOLE INSERT Ø 5.0
LSS-001IN**





Ortho Innovations.

PAKISTAN

CANNULATED LOCKING SCREW 5.0 Ø MM

CSL-301A	30MM
CSL-302A	35MM
CSL-303A	40MM
CSL-304A	45MM
CSL-305A	50MM
CSL-306A	55MM
CSL-307A	60MM
CSL-308A	65MM
CSL-309A	70MM
CSL-310A	75MM
CSL-311A	80MM
CSL-312A	85MM
CSL-313A	90MM
CSL-314A	95MM
CSL-315A	100MM
CSL-316A	105MM
CSL-317A	110MM
CSL-318A	115MM
CSL-319A	120MM
CSL-320A	125MM
CSL-321A	130MM
CSL-322A	135MM
CSL-323A	140MM
CSL-324A	145MM



IN TITANIUM CSL-301B TO CSL-324B

CANNULATED LOCKING SCREW 7.3Ø MM

CSL-301C	30MM
CSL-302C	35MM
CSL-303C	40MM
CSL-304C	45MM
CSL-305C	50MM
CSL-306C	55MM
CSL-307C	60MM
CSL-308C	65MM
CSL-309C	70MM
CSL-310C	75MM
CSL-311C	80MM
CSL-312C	85MM



OrthoInnovations.

PAKISTAN

<i>CSL-313C</i>	<i>90MM</i>
<i>CSL-314C</i>	<i>95MM</i>
<i>CSL-315C</i>	<i>100MM</i>
<i>CSL-316C</i>	<i>105MM</i>
<i>CSL-317C</i>	<i>110MM</i>
<i>CSL-318C</i>	<i>115MM</i>
<i>CSL-319C</i>	<i>120MM</i>
<i>CSL-320C</i>	<i>125MM</i>
<i>CSL-321C</i>	<i>130MM</i>
<i>CSL-322C</i>	<i>135MM</i>
<i>CSL-323C</i>	<i>140MM</i>
<i>CSL-324C</i>	<i>145MM</i>
<i>IN TITANIUM CSL-301D TO CSL-324D</i>	

