




Ortho Innovations.

PAKISTAN

ORTHON MINI INSTRUMENTS SET





Ortho Innovations.

PAKISTAN

JACOB CHUCK DRILL BITS

MIS-001	1.1MM \emptyset	45MM LENGTH
MIS-002	1.5MM \emptyset	70MM LENGTH
MIS-003	2.0MM \emptyset	85MM LENGTH
MIS-004	2.7MM \emptyset	85MMLENGTH



QUICK COUPLING DRILL BITS

MIS-005	1.1MM \emptyset	60MM LENGTH
MIS-006	1.5MM \emptyset	85MM LENGTH
MIS-007	1.5MM \emptyset	110MM LENGTH
MIS-008	2.0MM \emptyset	100MM LENGTH
MIS-009	2.0MM \emptyset	125MM LENGTH
MIS-010	2.7MM \emptyset	100MM LENGTH
MIS-011	2.7MM \emptyset	125MM LENGTH



QUICK COUPLINGTAP

MIS-012	1.5MM \emptyset	50MM LENGTH
MIS-013	2.0MM \emptyset	53MM LENGTH
MIS-014	2.0MM \emptyset	100MM LENGTH
MIS-015	2.7MM \emptyset	100MM LENGTH
MIS-016	2.7MM \emptyset	125MM LENGTH



MIS-017 QUICK COUPLING SCREW DRIVER SHAFT



MIS-018 HEXAGONAL SCREW DRIVER WITH GROOVE





Ortho Innovations.

PAKISTAN

MIS-019 HOLDING SLEEVE FOR MIS-O-018



MIS-020 MINI SCREW DRIVER FOR
1.5MM-2.0MM SCREWS



MIS-021 COUNTERSINK



MIS-022 MINI QUICK COUPLING HANDLE



MIS-023 QUICK COUPLING HANDLE



MIS-024 DOUBLE DRILL SLEEVE 1.1MM/1.5MM



MIS-025 DOUBLE DRILL SLEEVE 1.5MM/2.0MM





Ortho Innovations.

PAKISTAN

MIS-026 NEUTRAL & LOAD DRILL SLEEVE FOR 1.5MM
DRILL BITS



MIS-027 NEUTRAL & LOAD DRILL SLEEVE FOR 2.0MM
DRILL BITS



MIS-028 SHORT DRILL GUIDE FOR 2.7MM PLATES



MIS-029 LONG DRILL GUIDE FOR 2.7MM PLATES



MIS-030 MINI DEPTH GAUGE FOR
1.5MM-2.0MM SCREWS



MIS-031 SMALL DEPTH GAUGE FOR
2.7MM-4.0MM SCREWS





Ortho Innovations

PAKISTAN

MIS-032 SHARP HOOK



MIS-033 SELF HOLDING SCREW FORCEPS



MIS-034 BENDING IRON FOR SMALL PLATES



MIS-035 SMALL RETRACTOR 6MM WIDE
160MM LENGTH SHORT NARROW TIP



MIS-036 SMALL RETRACTOR 6MM WIDE
160MM LENGTH SHORT NARROW TIP



MIS-037 SMALL RETRACTOR 15MM



MIS-038 PERIOSTEAL ELEVATOR 3MM WIDE





Ortho Innovations.

PAKISTAN

MIS-039 MINI PLATEHOLDING FORCEPS



MIS-040 SMALL REDUCTION FORCEPS POINTED
130MM LENGTH



MIS-041 K-WIRES BENDING IRON 120MM LENGTH



MIS-042 MULTIPLE HOLE DRILL SLEEVE



MIS-043 WIRE BENDING PLIER 155MM LENGTH





ORTHON MINI IMPLANTS SET

2.0MM MINI DYNAMIC COMPRESSION PLATES

CAT NO	HOLES	LENGTH
--------	-------	--------

MCP-001	4	22MM
---------	---	------

MCP-002	5	27MM
---------	---	------

MCP-003	6	32MM
---------	---	------

MCP-004	7	37MM
---------	---	------

MCP-005	8	42MM
---------	---	------

IN TITANIUM FROM MCP-2001 TO MCP-2005



PROFILE

THICKNESS 1.3MM

WIDTH 5.1MM

HOLE SPACING 5.0MM

2.0MM MINI L & R PLATES

MCP-101	LEFT
---------	------

MCP-102	RIGHT
---------	-------

MCP-103	LEFT OBLIQUE
---------	--------------

MCP-104	RIGHT OBLIQUE
---------	---------------

MCP-105	T-SHAPED
---------	----------

IN TITANIUM MCP-3001 TO MCP-3005



PROFILE

THICKNESS 1.0MM

WIDTH 5.0MM

1.5MM MINI STRAIGHT PLATE

MCP-201	20	100MM
---------	----	-------

IN TITANIUM MCP-4001



2.0MM MINI ADAPTION PLATE

MCP-301	20	100MM
---------	----	-------

IN TITANIUM MCP-5001





1.5MM MINI T-PLATES

MCP-401	3/3
MCP-402	3/5
MCP-403	4/4
MCP-404	4/6

IN TITANIUM MCP-6001 TO MCP-6004



2.0MM MINI T-PLATES

MCP-501	3/3
MCP-502	3/5
MCP-503	4/4
MCP-504	4/6

IN TITANIUM MCP-7001 TO MCP-7005



2.7MM STRAIGHT RECONSTRUCTION PLATES

MCP-601	6	48MM
MCP-602	8	64MM
MCP-603	10	80MM
MCP-604	12	96MM
MCP-605	14	112MM
MCP-606	16	128MM

IN TITANIUM MCP-8001 TO MCP-8006



<i>PROFILE</i>	
<i>THICKNESS</i>	3.2MM
<i>WIDTH</i>	8.0MM
<i>HOLE SPACING</i>	8.0MM

2.7MM MINI DYNAMIC COMPRESSION PLATES

MCP-701	2	20MM
MCP-702	3	28MM
MCP-703	4	36MM
MCP-704	5	44MM
MCP-705	6	52MM
MCP-706	7	60MM
MCP-707	8	68MM
MCP-708	9	76MM
MCP-709	10	84MM
MCP-710	11	92MM
MCP-711	12	100MM

IN TITANIUM MCP-9001 TO MCP-9011



<i>PROFILE</i>	
<i>THICKNESS</i>	2.3MM
<i>WIDTH</i>	8.1MM
<i>HOLE SPACING</i>	8.0&12.0MM



1/4th TUBULAR PLATES

MCP-801	3	23MM
MCP-802	4	31MM
MCP-803	5	39MM
MCP-804	6	47MM
MCP-805	7	55MM
MCP-806	8	63MM

IN TITANIUM MCP-10001 TO MCP-10008



<i>PROFILE</i>	
THICKNESS	1.0MM
WIDTH	7.0MM
HOLE SPACING	8.0MM

2.7MM MINI L&T PLATES

MCP-901	LEFT
MCP-902	RIGHT
MCP-903	LEFT OBLIQUE
MCP-904	RIGHT OBLIQUE
MCP-905	T-SHAPED

IN TITANIUM MCP-10011 TO MCP-10015



<i>PROFILE</i>	
THICKNESS	1.2MM
WIDTH	7.0MM

K-WIRES

MCP-1001	.8MM	70MM
MCP-1002	1.0MM	150MM
MCP-1003	1.2MM	150MM
MCP-1004	1.4MM	150MM
MCP-1005	1.5MM	150MM
MCP-1006	1.6MM	150MM
MCP-1007	1.8MM	150MM
MCP-1008	2.0MM	150MM
MCP-1009	2.5MM	150MM
MCP-1010	3.0MM	150MM

IN TITANIUM MCP-10111 TO MCP-10121





Ortho Innovations.

PAKISTAN

SUTURE WIRES

MCP-1101	.1MM	10METERS
MCP-1102	.2MM	10METERS
MCP-1103	.3MM	10MWTERS
MCP-1104	.4MM	10MWTERS
MCP-1105	.5MM	10METERS
IN TITANIUM MCP-11111 TO MCP-11115		





Ortho Innovations.

PAKISTAN

ORTHON MINI SCREW SET

1.5MM CORTEX SCREWS

CAT NO	LENGTH
MSS-001	6MM
MSS-002	7MM
MSS-003	8MM
MSS-004	9MM
MSS-005	10MM
MSS-006	11MM
MSS-007	12MM
MSS-008	14MM
MSS-009	16MM
MSS-010	18MM
MSS-011	20MM



PROFILE	
HEAD DIAMETER	3.0MM
HEXAGONAL SOCKET WIDTH	1.5MM
THREAD DIAMETER	1.5MM
CORE DIAMETER	1.0MM
DRILL BIT FOR THREADED HOLES	1.1MM
DRILL BIT FOR GLIDING HOLES	1.5MM
TAP	1.5MM

IN TITANIUM MSS-401 TO MSS-411

2.0MM CORTEX SCREWS

MSS-102	6MM
MSS-103	7MM
MSS-104	8MM
MSS-105	9MM
MSS-106	10MM
MSS-107	11MM
MSS-108	12MM
MSS-109	14MM
MSS-110	16MM
MSS-111	18MM
MSS-112	20MM
MSS-113	22MM
MSS-114	24MM
MSS-115	26MM
MSS-116	28MM
MSS-117	30MM
MSS-118	32MM
MSS-119	34MM



PROFILE	
HEAD DIAMETER	4.0MM
HEXAGONAL SOCKET WIDTH	1.5MM
THREAD DIAMETER	2.0MM
CORE DIAMETER	1.3MM
DRILL BIT FOR THREADED HOLES	1.5MM
DRILL BIT FOR GLIDING HOLES	2.0MM
TAP	2.0MM



MSS-120 36MM
MSS-121 38MM
MSS-122 40MM
IN TITANIUM MSS-501 TO MSS-522

2.7MM CORTEX SCREWS

MSS-201 6MM
MSS-202 7MM
MSS-203 8MM
MSS-204 9MM
MSS-205 10MM
MSS-206 11MM
MSS-207 12MM
MSS-208 14MM
MSS-209 16MM
MSS-210 18MM
MSS-211 20MM
MSS-212 22MM
MSS-213 24MM
MSS-214 26MM
MSS-215 28MM
MSS-216 30MM
MSS-217 32MM
MSS-218 34MM
MSS-219 36MM
MSS-220 38MM
MSS-221 40MM
IN TITANIUM MSS-601 TO MSS-621



PROFILE	
HEAD DIAMETER	5.0MM
HEXAGONAL SOCKET WIDTH	2.5MM
THREAD DIAMETER	2.7MM
CORE DIAMETER	1.9MM
DRILL BIT FOR THREADED HOLES	2.0MM
DRILL BIT FOR GLIDING HOLES	2.7MM
TAP	2.7MM

SMALL WASHER 7.0MM DIAMETER

MSS-701

