




Ortho Innovations.

PAKISTAN

PELVIC REDUCTION INSTRUMENTS SET





Ortho Innovations

PAKISTAN

PIS-001

JACOB CHUCK



PIS-002

KEYLESS CHUCK WITH T-HANDLE



PIS-003

BONE HOOK LARGE
200MM LONG



PIS-004

SHARP BONE HOOK



PIS-005

PELVIC RETRACTOR BLUNT
270MM LONG





Ortho Innovations.

PAKISTAN

PIS-006

BALL SPIKE WITH BALL TIP
STRAIGHT LENGTH 300MM



PIS-007

DISC, SPIKED, FOR PELVIC
REDUCTION FORCEP



PIS-008

PELVIC REDUCTION FORCEPS
SMALL, LENGTH 200MM



PIS-009

PELVIC REDUCTION FORCEPS
MEDIUM, LENGTH 250MM





Ortho Innovations.

PAKISTAN


PIS-010

PELVIC REDUCTION FORCEPS
LARGE, LENGTH 300MM



PIS-011

PELVIC REDUCTION FORCEPSSMALL,
OBLIQUE, WITH POINTED TIP, LENGTH 190MM



PIS-012

PELVIC REDUCTION FORCEPS LARGE
OBLIQUE, WITH POINTED BALL TIPS, LENGTH 240MM





Ortho Innovations.

PAKISTAN

PIS-013

PELVIC REDUCTION FORCEPS, LONG
WITH POINTED BALL TIPS, LENGTH 400MM



PIS-014

PELVIC REDUCTION FORCEPS, LONG
WITH 3 POINTED BALL TIPS, LENGTH 400MM



PIS-015

PELVIC REDUCTION FORCEPS WITH ANGLED
POINTED BALL TIPS, LENGTH 250MM





Ortho Innovations.

PAKISTAN

PIS-016

LARGE PINTED REDUCTION FORCEPS
LENGTH, 200MM



PIS-017

BENDING IRON IN PAIRS
LENGTH, 190MM



Ortho Innovations.


PAKISTAN

IMPLANTS (BOTH STAINLESS STEEL AND TITANIUM)



3.5MM RECONSTRUCTION PLATES STRAIGHT

CAT NO	HOLES	LENGTH
PIM-001	3	34MM
PIM-002	4	46MM
PIM-003	5	70MM
PIM-004	6	84MM
PIM-005	7	99MM
PIM-006	8	113MM
PIM-007	9	128MM
PIM-008	10	142MM
PIM-009	12	171MM
PIM-010	14	200MM
PIM-011	16	228MM
PIM-012	18	257MM
PIM-013	20	286MM
PIM-014	22	315MM



IN TITANIUM PIM101 TO PIM-113

LCP RECONSTRUCTION PLATE 3.5MM COMBINED HOLES

PIM-001A TO PIM-014 AND PIM/101A TO PIM-113A IN TITANIUM



3.5MM RECONSTRUCTION PLATES CURVED

PIM-201	6	70MM
PIM-202	8	94MM
PIM-203	10	118MM
PIM-204	12	142MM
PIM-205	14	166MM
PIM-206	16	190MM
PIM-207	18	214MM

IN TITANIUM PIM-301TO PIM-307



<i>PROFILE</i>	
<i>THICKNESS</i>	<i>2.8MM</i>
<i>WIDTH</i>	<i>10.0MM</i>
<i>HOLE SPACING</i>	<i>12.0MM</i>

4.5MM RECONSTRUCTION PLATES STRAIGHT

PIM-401	3	45MM
PIM-402	4	61MM
PIM-403	5	77MM
PIM-404	6	93MM
PIM-405	7	109MM
PIM-406	8	125MM
PIM-407	9	141MM
PIM-408	10	157MM
PIM-409	12	189MM
PIM-410	14	221MM
PIM-411	16	253MM

IN TITANIUM PIM-501TO PIM-511



<i>PROFILE</i>	
<i>THICKNESS</i>	<i>2.8MM</i>
<i>WIDTH</i>	<i>12.0MM</i>
<i>HOLE SPACING</i>	<i>16.0MM</i>



Ortho Innovations.

PAKISTAN

4.5MM RECONSTRUCTION PLATES CURVED

PIM-601	6	93MM
PIM-602	8	125MM
PIM-603	10	157MM
PIM-604	12	189MM
PIM-605	14	221MM

IN TITANIUM PIM-701 TO PIM-705



PROFILE	
THICKNESS	2.8MM
WIDTH	12.0MM
HOLE SPACING	16.0MM

SCHANZ SCREW 6.0MM/190MM

PIS-801



PIS-910 NARROW DCP, 2 HOLES
LENGTH, 39MM



LSS-101	25MM
LSS-102	30MM
LSS-103	35MM
LSS-104	40MM
LSS-105	45MM
LSS-106	50MM
LSS-107	55MM
LSS-108	60MM
LSS-109	65MM
LSS-110	70MM
LSS-111	75MM
LSS-112	80MM
LSS-113	85MM
LSS-114	90MM





Ortho Innovations.

PAKISTAN

LSS-115	95MM
LSS-116	100MM
LSS-117	105MM
LSS-118	110MM
LSS-119	115MM
LSS-120	120MM

6.5MM CANCELLOUS SCREW 32MM THREADS

LSS-122	25MM
LSS-123	30MM
LSS-124	35MM
LSS-125	40MM
LSS-126	45MM
LSS-127	50MM
LSS-128	55MM
LSS-129	60MM
LSS-130	65MM
LSS-131	70MM
LSS-132	75MM
LSS-133	80MM
LSS-134	85MM
LSS-135	90MM
LSS-136	95MM
LSS-137	100MM
LSS-138	105MM
LSS-139	110MM
LSS-140	115MM
LSS-141	120MM

PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>SHAFT DIAMETER</i>	<i>4.5MM</i>
<i>THREAD DIAMETER</i>	<i>6.5MM</i>
<i>CORE DIAMETER</i>	<i>3.0MM</i>
<i>THREAD LENGTH</i>	<i>16MM</i>
<i>DRILL BIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>6.5MM</i>



PROFILE

<i>HEAD DIAMETER</i>	<i>8.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>3.5MM</i>
<i>SHAFT DIAMETER</i>	<i>4.5MM</i>
<i>THREAD DIAMETER</i>	<i>6.5MM</i>
<i>CORE DIAMETER</i>	<i>3.0MM</i>
<i>THREAD LENGTH</i>	<i>32MM</i>
<i>DRILLBIT</i>	<i>3.2MM</i>
<i>TAP</i>	<i>6.5MM</i>



6.5MM CANCELLOUS SCREW FULL THREAD

LSS-142	25MM
LSS-143	30MM
LSS-144	35MM
LSS-145	40MM
LSS-146	45MM
LSS-147	50MM
LSS-148	55MM
LSS-149	60MM
LSS-150	65MM
LSS-151	70MM
LSS-152	75MM
LSS-153	80MM
LSS-154	85MM
LSS-155	90MM
LSS-156	95MM



PROFILE	
HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	6.5MM
DRILL BIT	3.2MM
TAP	6.5MM

4.5MM CORTEX SCREWS SELF TAPPING

LSS-056	14MM
LSS-057	16MM
LSS-058	18MM
LSS-059	20MM
LSS-060	22MM
LSS-061	24MM
LSS-062	26MM
LSS-063	28MM
LSS-064	30MM
LSS-065	32MM
LSS-066	34MM
LSS-067	36MM
LSS-068	38MM
LSS-069	40MM
LSS-070	42MM
LSS-071	44MM
LSS-072	46MM
LSS-073	48MM
LSS-074	50MM





LSS-075	52MM
LSS-076	54MM
LSS-077	56MM
LSS-078	58MM
LSS-079	60MM
LSS-080	62MM
LSS-081	64MM
LSS-082	66MM
LSS-083	68MM
LSS-084	70MM
LSS-085	72MM
LSS-086	74MM
LSS-087	76MM
LSS-088	78MM
LSS-089	80MM
LSS-090	82MM
LSS-091	84MM
LSS-092	86MM
LSS-093	88MM
LSS-094	90MM
LSS-095	95MM
LSS-096	100MM
LSS-097	105MM
LSS-098	110MM
LSS-099	115MM
LSS-100	120MM
IN TITANIUM LSS-156 TO LSS-200	

PROFILE	
HEAD DIAMETER	8.0MM
HEXAGONAL SOCKET WIDTH	3.5MM
THREAD DIAMETER	4.5MM
CORE DIAMETER	3.0MM
DRILL BIT	3.2MM

3.5MM CORETX SCREW SELF TAPPING

SSS-026	10MM
SSS-027	12MM
SSS-028	14MM
SSS-029	16MM
SSS-030	18MM
SSS-031	20MM
SSS-032	22MM
SSS-033	24MM
SSS-034	26MM
SSS-036	28MM
SSS-037	30MM





SSS-038	32MM
SSS-039	34MM
SSS-040	36MM
SSS-041	38MM
SSS-042	40MM
SSS-043	42MM
SSS-044	44MM
SSS-045	46MM
SSS-046	48MM
SSS-047	50MM
SSS-048	55MM
SSS-049	60MM
SSS-050	65MM
SSS-051	70MM

IN TITANIUM SSS-106 TO SSS-151

(FOR LOCKING HEAD SCREWS SEE SECTION ON SCREWS)

PROFILE
HEAD DIAMETER 6.0MM
HEXAGONAL SOCKET WIDTH 2.5MM
THREAD DIAMETER 3.5MM
CORE DIAMETER 2.4MM
DRILL BIT FOR GLIDING HOLES
3.5MM
DRILL BIT FOR THREADED HOLES

PIM-606 SYMPHYSEAL PLATE 3.5 WITH COMBI HOLES
6 HOLES 78MM LENGTH



PIM-607 SYMPHYSEAL PLATE 3.5 WITH COMBI HOLES
AND 2 DCP 6 HOLES 78MM LENGTH



