




Orthon Innovations.

PAKISTAN

ORTHON SMALL INSTRUMENTS SET



DRILL BIT

SIS-001	2.5MM ϕ	95MM LENGTH
SIS-002	3.5MM ϕ	95MM LENGTH
SIS-003	3.5MM ϕ	180MM LENGTH



QUICK COUPLING DRILL BIT

SIS-004	2.5MM ϕ	110MM LENGTH
SIS-005	2.5MM ϕ	180MM LENGTH
SIS-006	3.5MM ϕ	100MM LENGTH
SIS-007	3.5MM ϕ	195MM LENGTH



QUICK COUPLING TAP

SIS-008	3.5MM ϕ	110MM LENGTH
SIS-009	3.5MM ϕ	175MM LENGTH
SIS-010	4.0MM ϕ	110MM LENGTH





Ortho Innovations.

PAKISTAN

SIS-011

SMALL HEXAGONAL SCREW DRIVER SHAFT



SIS-012

SMALL HEXAGONAL SCREW DRIVER



SIS-013

SMALL HOLDING SLEEVE



SIS-014

SMALL COUNTERSINK



SIS-015

QUICK COUPLING HANDLE





Ortho Innovations.

PAKISTAN

SIS-016

QUICK COUPLING T-HANDLE



SIS-017

DOUBLE DRILL SLEEVE 3.5/2.5MM



SIS-018

INSERT DRILL SLEEVE 3.5/2.5MM



SIS-O-019

SMALL DEPTH GAUGE



SIS-020

SHARP HOOK



SIS-021

SELF HOLDING SCREW FORCEPS





Ortho Innovations.

PAKISTAN

SIS-022

BENDING PLIER



SIS-023

BENDING IRON 4.5/2.5MM



SIS-024

BENDING IRON 2.5/4.5MM



SIS-025

DRILL SLEEVE (NEUTRAL&LOAD)
DRILL BIT 2.5MM ϕ



SIS-026

DRILL SLEEVE (NEUTRAL & LOAD)
FOR SMALL LIMITED CONTACT PLATES



SIS-027

UNIVERSAL DRILL SLEEVE 2.5/3.5MM
FOR SMALL LIMITED CONTACT PLATES



SIS-028

BENDING TEMPLATE 87MM LENGTH
BENDING TEMPLATE 114MM LENGTH





Ortho Innovations.

PAKISTAN

SIS-029

SMALL RETRACTER 6MM WIDE SHORT
NARROW TIP 160 MM LENGTHS



SIS-030

SMALL RETRACTER 8MM WIDE
SHORT NARROW TIP 160MM LENGTH



SIS-031

SMALL RETRACTER 15MM WIDE



SIS-032

BONE HOLDING FORCEPS
SELF CENTRING SIZE 0-190MM LENGTH



SIS-033

SMALL REDUCTION FORCEPS WITH POINTS
130MMLLENGTH





Ortho Innovations.

PAKISTAN

SIS-034

SMALL REDUCTION FORCEPS
SERRATED JAWS 140MM LENGTH



SIS-035

PERIOSTEAL ELEVATOR 6MM WIDE
SLIGHTLY CURVED BLADE, STRAIGHT EDGE



SIS-036

PERIOSTEAL ELEVATOR 6MM WIDE
SLIGHTLY CURVED BLADE, ROUND EDGE



SIS-037

K-WIRE BENDING IRON FOR 1.2TO2.5MM
LENGTH 120MM



SIS-038

TRIPPLE DRILL DUIDE FOR K-WIRES



SIS-039


WIRE BENDING PLIER 155MM LENGTH





1/3rd TUBULAR PLATES

CAT NO	HOLES	LENGTH
SCP-001	2	25MM
SCP-002	3	37MM
SCP-003	4	49MM
SCP-004	5	61MM
SCP-005	6	73MM
SCP-006	7	85MM
SCP-007	8	97MM
SCP-008	9	109MM
SCP-009	10	121MM
SCP-010	12	145MM



PFOFILE	
THICKNESS	1.0MM
WIDTH	9.0MM
HOLE SPACING	12.0&16.0MM

IN TITANIUM SCP-101 TO 110

SMALL DYNAMIC COMPRESSION PLATES

SCP-201	2	25MM
SCP-202	3	37MM
SCP-203	4	49MM
SCP-204	5	61MM
SCP-205	6	73MM
SCP-206	7	85MM
SCP-207	8	97MM
SCP-208	9	109MM
SCP-209	10	121MM
SCP-210	11	134MM
SCP-211	12	145MM



PROFILE	
THICKNESS	3.2MM
WIDTH	10.0MM
HOLE SPACING	12.0716.0MM

IN TITANIUM SCP-301 TO SCP-311

SMALL LIMITED CONTACT DYNAMIC COMPRESSION PLATES

SLCP-001	4	51MM
SLCP-002	5	64MM
SLCP-003	6	77MM
SLCP-004	7	90MM
SLCP-005	8	103MM
SLCP-006	9	116MM
SLCP-007	10	129MM
SLCP-008	11	142MM
SLCP-009	12	155MM



PROFILE	
THICKNESS	4.1MM
WIDTH	11.0MM
HOLE SPACING	13.0MM


IN TITANIUM SLCP-101 TO SLCP-109



SMALL T – PLATES

SCP-401	3/3	39MM
SCP-402	4/4	56MM
SCP-403	5/3	67MM
SCP-404	6/4	78MM

IN TITANIUM SCP-501 TO SCP-504



<i>PROFILE</i>	
THICKNESS	1.3MM
WIDTH	10.0MM

SMALL T – PLATES, OBLIQUE


SCP-601	3/3	52MM
SCP-602	4/3	63MM
SCP-603	5/3	74MM

IN TITANIUM SCP-701 TO SCP-703

LCP-T-PLATE 3.5MM OBLIQUE-ANGLED

SCP-601A TO SCP-603A FOR RIGHT
 SCP-601B TO SCP-603B FOR LEFT
 SCP-601C TO SCP-603C FOR R.TITANIUM
 SCP-601D TO SCP-603D FOR L.TATINIUM

CLOVERLEAF PALTES



<i>PROFILE</i>	
THICKNESS	1.6MM
WIDTH	10.0MM



<i>PROFILE</i>	
THICKNESS	2.1MM
WIDTH	15.2MM
HEAD THICKNESS	1.2MM
SHAFT WIDTH	15.3MM

SCP-801	3	88MM
SCP-802	4	104MM
SCP-803	5	120MM
SCP-804	6	136MM

IN TITANIUM SCP-901 TO SCP-904

LCP CLOVER LEAF

SCP-801A TO SCP-804A
 SCP-801B TO SCP-804B IN TITANIUM

Y – PLATES

SCP-1001	4
SCP-1002	6
SCP-1003	8



<i>PROFILE</i>	
THICKNESS	2.0MM

IN TITANIUM SCP-2001 TO SCP-2003



Ortho Innovations.

PAKISTAN

CALCANEAL PLATES

SCP-3001 60MM

SCP-3002 70MM

IN TITANIUM SCP-4001 TO SCP-4002



PROFILE
THICKNESS 1.2MM

CALCANEAL PLATES – (Y – SHAPE)

SCP-5001 80MM

IN TITANIUM SCP-6001



PROFILE
THICKNESS 1.2MM

KIRSCHNER WIRES WITH SINGLE TROCER TIP

SCP-7001 0.8MM 70MM

SCP-7002 1.0MM 150MM

SCP-7003 1.2MM 150MM

SCP-7004 1.4MM 150MM

SCP-7005 1.5MM 150MM

SCP-7006 1.6MM 150MM

SCP-7007 1.8MM 150MM

SCP-7008 2.0MM 150MM

SCP-7009 2.5MM 150MM

SCP-7010 3.0MM 150MM





LCP HOOK PLATE 3.5MM 3 HOLES 62MM LENGTH

- SCP-8001 SS
- SCP-9001 TITANIUM



LCP METAPHYSEAL PLATE 3.5MM

CAT NO	HOLES	LENGTH
SCP-10001	6	86
SCP-10002	7	99
SCP-10003	8	112
SCP-10004	9	125
SCP-10005	10	138
SCP-10006	11	151
SCP-10007	12	164
SCP-10008	14	190
SCP-10009	16	216
SCP-10010	18	242

IN TITANIUM SCP-20001 TO SCP-20010



LCP 3.5 MM

CAT NO	HOLES
SCP-30001	4
SCP-30002	5
SCP-30003	6
SCP-30004	7
SCP-30005	8
SCP-30006	9
SCP-30007	10
SCP-30008	11

IN TITANIUM SCP-40001 TO 40008



PROFILE
THICKNESS 3.3MM
HEIGHT 4.0MM
WIDTH 11.0MM
HOLE SPACING 13.0MM



LCP PILON PLATE 2.7/3.5 CRUCIFORM

CAT NO	HOLES SHAFT
SCP-50001	7
SCP-50002	9
IN TITANIUM SCP-60001 TO SCP-60002	



LCP-T PLATE 3.5MM RIGHT ANGLED

CAT NO	HOLES HEAD	HOLES SHAFT	LENGTH
SCP-70001	3	3	50
SCP-70002	4	4	56
SCP-70003	3	5	67
SCP-70004	4	6	78
IN TITANIUM SCP-80001 TO SCP-80004			



CALCANEAL LOCKING PLATE 3.5MM

CAT NO	DIRECTION	SIZE	LENGTH
SCP-90001	RIGHT	EXTRA SMALL	64MM
SCP-90002	LEFT	EXTRA SMALL	64MM
SCP-90003	RIGHT	SMALL	69MM
SCP-90004	LEFT	SMALL	69MM
SCP-90005	RIGHT	LARGE	76MM
SCP-90006	LEFT	LARGE	76MM
SCP-90007	RIGHT	EXT.LARGE	81MM
SCP-90008	LEFT	EXT.LARGE	81MM



IN TITANIUM SCP-100001 TO SCP-100008



Ortho Innovations.

PAKISTAN

LCP SPACER Ø 3.5MM, LENGTH 2.0MM

ALLOWS THE 2.0MM PLATE TO BE LIFTED FROM BONE SURFACE
TO PRESERVE VASCULARITY.

SCP-1000001 AND IN TITANIUM SCP-100001



WASHER Ø 3.5/2.7MM

ALLOWS INSERTION OF CORTEX SCREW Ø2.7MM IN THE DC PART OF THE LCP 3.5

SCP-1000001A AND IN TITANIUM SCP-1000001B

