



*Ortho Innovations.*

PAKISTAN

## ORTHON SMALL SCREW SET



### 3.5MM CORTEX SCREWS IN STAINLESS STEEL AND TITANIUM

CAT NO	SIZE
SSS-001	10MM
SSS-002	12MM
SSS-003	14MM
SSS-004	16MM
SSS-005	18MM
SSS-006	20MM
SSS-007	22MM
SSS-008	24MM
SSS-009	26MM
SSS-010	28MM
SSS-011	30MM
SSS-012	32MM
SSS-013	34MM
SSS-014	36MM
SSS-015	38MM
SSS-016	40MM
SSS-017	42MM
SSS-018	44MM
SSS-019	46MM



PROFILE  
HEAD DIAMETER 6.0MM  
HEXAGONAL SOCKET WIDTH 2.5MM  
THREAD DIAMETER 3.5MM  
CORE DIAMETER 2.4MM  
DRILL BIT FOR GLIDING HOLES  
3.5MM  
DRILL BIT FOR THREADED HOLES



*Ortho Innovations*

PAKISTAN

SSS-020	48MM
SSS-021	50MM
SSS-022	55MM
SSS-023	60MM
SSS-024	65MM
SSS-025	70MM

**IN TITANIUM SSS-101 TO SSS-125**

**3.5MM CORETX SCREW SELF TAPPING IN STAINLESS STEEL AND TITANIUM**

SSS-201	10MM
SSS-202	12MM
SSS-203	14MM
SSS-204	16MM
SSS-205	18MM
SSS-206	20MM
SSS-207	22MM
SSS-208	24MM
SSS-209	26MM
SSS-210	28MM
SSS-211	30MM
SSS-212	32MM
SSS-213	34MM
SSS-214	36MM
SSS-215	38MM
SSS-216	40MM
SSS-217	42MM
SSS-218	44MM
SSS-219	46MM
SSS-220	48MM
SSS-221	50MM
SSS-222	55MM
SSS-223	60MM
SSS-224	65MM
SSS-225	70MM

**IN TITANIUM SSS-301 TO SSS-325**



*PROFILE*

<i>HEAD DIAMETER</i>	<i>6.0MM</i>
<i>HEXAGONAL SOCKET WIDTH</i>	<i>2.5MM</i>
<i>THREAD DIAMETER</i>	<i>3.5MM</i>
<i>CORE DIAMETER</i>	<i>2.4MM</i>

*DRILL BIT FOR GLIDING HOLES*  
*3.5MM*

*DRILL BIT FOR THREADED HOLES*



**4.0MM CANCELLOUS SCREW SHORT THREADED IN STAINLESS STEEL AND TITANIUM**

SSS-401	10MM
SSS-402	12MM
SSS-403	14MM
SSS-404	16MM
SSS-405	18MM
SSS-406	20MM
SSS-407	22MM
SSS-408	24MM
SSS-409	26MM
SSS-410	28MM
SSS-411	30MM
SSS-412	32MM
SSS-413	34MM
SSS-414	36MM
SSS-415	38MM
SSS-416	40MM
SSS-417	42MM
SSS-418	45MM
SSS-419	50MM
SSS-420	55MM
SSS-421	60MM
SSS-422	65MM
SSS-423	70MM



PROFILE	
HEAD DIAMETER	6.0MM
HEXAGONAL SOCKET WIDTH	2.5MM
THREAD DIAMETER	4.0MM
CORE DIAMETER	1.9MM
DRILL BIT FOR THREADED HOLES	2.0MM
TAP	4.0MM

**IN TITANIUM SSS-501 TO SSS-523**

**4.0MM CANCELLOUS SCREW FULLY THREADED**

SSS-601	10MM
SSS-602	12MM
SSS-603	14MM
SSS-604	16MM
SSS-605	18MM
SSS-606	20MM
SSS-607	22MM
SSS-608	24MM
SSS-609	26MM
SSS-610	28MM





SSS-611	30MM
SSS-612	32MM
SSS-613	34MM
SSS-614	36MM
SSS-615	38MM
SSS-616	40MM
SSS-617	42MM
SSS-618	45MM
SSS-619	50MM
SSS-620	55MM
SSS-621	60MM
SSS-622	65MM
SSS-623	70MM

**IN TITANIUM SSS-701 TO SSS-723**

PROFILE	
HEAD DIAMETER	6.0MM
HEXAGONAL SOCKET WIDTH	2.5MM
THREAD DIAMETER	4.0MM
CORE DIAMETER	1.9MM
DRILL BIT FOR THREADED HOLES	2.0MM
TAP	4.0MM


### WASHER SMALL 7.0MM SIZE

SSS-801  
**IN TITANIUM SSS-901**



### LOCKING SCREW Ø 2.7MM (HEAD LCP 2.4)

CAT NO	SIZE
SSS-001A	10MM
SSS-002A	12MM
SSS-003A	14MM
SSS-004A	16MM
SSS-005A	18MM
SSS-006A	20MM
SSS-007A	22MM
SSS-008A	24MM
SSS-009A	26MM
SSS-010A	28MM
SSS-011A	30MM



PROFILE	
THREAD DIAMETER	2.7MM
DRILL BIT FOR THREADED HOLE	2.0MM
CORE DIAMETER	2.1MM
HEAD DIAMETER	3.5MM



*Ortho Innovations.*

PAKISTAN


SSS-012A	32MM
SSS-013A	34MM
SSS-014A	36MM
SSS-015A	38MM
SSS-016A	40MM
SSS-017A	42MM
SSS-018A	44MM
SSS-019A	46MM
SSS-020A	48MM
SSS-021A	50MM
SSS-022A	55MM
SSS-023A	60MM
SSS-024A	65MM
SSS-025A	70MM

**IN TITANIUM SSS-001B TO SSS-025B**

**LOCKING SCREW Ø 2.7MM (HEAD LCP3.5) SELF TAPPING**

SSS-006C	20MM
SSS-007C	22MM
SSS-008C	24MM
SSS-009C	26MM
SSS-010C	28MM
SSS-011C	30MM
SSS-012C	32MM
SSS-013C	34MM
SSS-014C	36MM
SSS-015C	38MM
SSS-016C	40MM
SSS-017C	42MM
SSS-018C	44MM
SSS-019C	46MM
SSS-020C	48MM
SSS-021C	50MM
SSS-022C	55MM
SSS-023C	60MM


**IN TITANIUM SSS-006D TO SSS.023D**





**LOCKING SCREW Ø 3.5MM SELF TAPPING**

SSS-201A	10MM
SSS-202A	12MM
SSS-203A	14MM
SSS-204A	16MM
SSS-205A	18MM
SSS-206A	20MM
SSS-207A	22MM
SSS-208A	24MM
SSS-209A	26MM
SSS-210A	28MM
SSS-211A	30MM
SSS-212A	32MM
SSS-213A	34MM
SSS-214A	36MM
SSS-215A	38MM
SSS-216A	40MM
SSS-217A	42MM
SSS-218A	44MM
SSS-219A	46MM
SSS-220A	48MM
SSS-221A	50MM
SSS-222A	55MM
SSS-223A	60MM
SSS-224A	65MM
SSS-225A	70MM
SSS-226A	75MM
SSS-227A	80MM
SSS-228A	85MM
SSS-229A	90MM



<i>PROFILE</i>	
<i>HEAD DIAMETER</i>	<i>5.0MM</i>
<i>THREAD DIAMETER</i>	<i>3.5MM</i>
<i>CORE DIAMETER</i>	<i>2.9MM</i>
<i>DRILL BIT FOR THREADED HOLES</i>	<i>2.8MM</i>
<i>TAP</i>	<i>3.5MM</i>

**IN TITANIUM SSS-201B TO SSS-229B**



*Ortho Innovations.*

PAKISTAN

**LOCKING SCREW Ø 3.5MM SELF DRILLING**

SSS-201B	10MM
SSS-202B	12MM
SSS-203B	14MM
SSS-204B	16MM
SSS-205B	18MM
SSS-206B	20MM
SSS-207B	22MM
SSS-208B	24MM
SSS-209B	26MM
SSS-210B	28MM
SSS-211B	30MM
SSS-212B	32MM
SSS-213B	34MM
SSS-214B	36MM
SSS-215B	38MM
SSS-216B	40MM
SSS-217B	42MM
SSS-218B	44MM
SSS-219B	46MM
SSS-220B	48MM
SSS-221B	50MM
SSS-222B	55MM
SSS-223B	60MM
SSS-224B	65MM
SSS-225B	70MM
SSS-226B	75MM
SSS-227B	80MM
SSS-228B	85MM
SSS-229B	90MM



<i>PROFILE</i>	
<i>HEAD DIAMETER</i>	<i>5.0MM</i>
<i>THREAD DIAMETER</i>	<i>3.5MM</i>
<i>CORE DIAMETER</i>	<i>2.9MM</i>
<i>DRILL BIT FOR THREADED HOLES</i>	<i>2.8MM</i>
<i>TAP</i>	<i>3.5MM</i>

**IN TITANIUM SSS-201C TO SSS-229C**