



*Ortho Innovations*

PAKISTAN

# ORTHON USS FRACTURE SYSTEM





**TRANSPEDICULAR SCHANZ SCREWS DUAL CORE**



THREAD LENGTH 35-55MM

Ø 5.0MM (USSF-001 TO USSF-005)

Ø 6.2MM (USSF-101 TO USSF-105)

**FRACTURE CLAMP FOR ROD**

Ø 6.0MM LOW PROFILE

USSF-201



**FRACTURE CLAMP FOR ROD**

Ø 6.0MM LOW PROFILE CRANIAL END

USSF-301



**ROD Ø 6.0MM HARD**

LENGTH 50MM TO 150MM

USSF-401 TO USSF-405



**CROSS LINK CLAMP FOR RODS**

Ø 6.0MM PREASSEMBLED

USSF-501



**ROD 3.5MM Ø FOR CROSS LINK STABILIZER**

LENGTH 40.0MM TO 80.0MM (USSF-601 TO USSF-605)





*Ortho Innovations*

PAKISTAN

**FIXATION RING FOR RODS**

Ø 6.0MM (USSF-701)



**TRANSPEDICULAR SCREWS WITH DUAL CORE AND DOUBLE THREAD, FOR REDUCTION OF SPONDYLITIS.**

THREAD LENGTH 40MM-50MM

Ø6.0MM (USSF-801 TO USSF-803)

Ø 7.0MM (USSF-901 TO USSF-903)

**ALL IMPLANTS ARE AVAILABLE IN BOTH TITANIUM AND STAINLESS STEEL.**

**FOR TITANIUM USSFT MAY BE USED INSTEAD OF USSF.**



**PRINCIPLE OF FRACTURE CLAMPS.**

**Controlled reduction due to the free angular play of  $\pm 18^\circ$   
(Fracture clamp for cranial end:  $+9^\circ-18^\circ$ )**

